**P65.3**

# Drehbare Aluminiumscheibe

## Rotatable Aluminium Disc

### Disque tournant en aluminium

#### 1. Beschreibung (Bild 1)

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#### 1. Description (Fig. 1)

The rotatable aluminium disc is used for observing the effect of eddy currents and, by making use of this characteristic, to demonstrate an A.C. meter.

The apparatus consists of an aluminium disc 150 mm in diameter fitted in a holder fork, and a  $\sqcup$ -shaped iron holder 30 mm wide. The holder fork runs into a supporting bar to be clamped to stand equipment.

#### 2. Arrangement to observe the effect of eddy currents

For this purpose the iron holder is not required and the disc is fitted to stand equipment so that it can be freely rotated between the plane surfaces of the pole-pieces (560 31) mounted on the U-core (562 11), namely in the magnetic field between these pole-pieces (see Directions for Use 560 31).

For further details on this experiment see: „Der zerlegbare Transformator“ (562 40) Experiment No. 31. (Only available in German.)

#### Notes

- The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
- The data DC... refer to the experiment descriptions published as a collection of leaflets in the „LEYBOLD PHYSIKALISCHE HANDBLÄTTER“ (PHYSICS LEAFLETS).
- The data and illustrations are not binding in every detail for the design of the apparatus. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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### 3. Arrangement to demonstrate an A. C. meter (Fig. 2)

For this purpose the rotatable aluminium disc is fitted to stand equipment in such a manner that it can rotate freely before a horizontal U-core (562 11) with two identical coils each with 250 turns (562 13), and between the iron holder with inserted transformer yoke (562 11 A) and the fitted coil with 10,000 turns (562 16).

For further details on this experiment see: "Der zerlegbare Transformator" (562 40) Experiment No. 74 to 79. (Only available in German.)

### 4. Other experiments

In a similar manner the rotatable aluminium disc is used in conjunction with the yoke (562 11 A) and the various types of coils (562 13 to 562 18) of the demonstration transformer to carry out experiments on "Rotation of a metal disc in an alternating magnetic field".

The disc can also be used as a diffraction disc for centimetre waves (DC 538.566; f).

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