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
140PC Series

Absolute, Differential, Gage, Vacuum Gage Sensors/Amplified

114

FEATURES

- FCB terminals on opposite side from the ports
- Optional color coded leadwires, 12 in., 24 gauge
- Fully signal conditioned

140PC Series PERFORMANCE CHARACTERISTICS at 8.0 ± 0.01 VDC Excitation, 25°C

	Min.	Typ.	Max.	Units
Excitation *	7.00	8.00	16.0	VDC
Supply Current	---	8.00	20.0	mA
Current Sourcing Output	---	---	10	mA
Null Offset (141/142PC)	0.95	1.00	1.05	V
Null Offset (143PC)**	3.45	3.50	3.55	V
Output at Full Pressure	5.90	6.00	6.10	V
Full Scale Output, F.S.O.† (141/142PC)	4.95	5.00	5.05	V
Full Scale Output, F.S.O.† (143PC)**	---	5.00	---	V
Ratiometricity Error				
7 to 8 V or 8 to 9 V	---	± 0.50	---	%F.S.O.
9 to 12 V	---	± 2.00	---	
Stability over One Year	---	± 0.50	---	%F.S.O.
Response Time	---	---	1.00	msec
Common Mode Pressure***	---	---	60	psi
Weight	---	5	---	grams
Short Circuit Protection	Output may be shorted indefinitely to ground			
Output Ripple	None, DC device			
Ground Reference	Supply and output are common			

* 8.0 VDC excitation is recommended with 1 psi unit.

** Positive and negative pressure measurement.

*** Higher common mode pressures possible if sensor is not used over entire operating temperature range.

† F.S.O. is defined as the algebraic difference between end points. Please note, actual output is 1 V to 6 V (at 8.00 ± 0.01 VDC). F.S.O. is then 5V.

ENVIRONMENTAL SPECIFICATIONS

Operating Temperature	-40° to +85°C (-40° to +185°F)
Storage Temperature	-65° to +125°C (-67° to +257°F)
Compensated Temperature	--- 18° to +63°C (0° to +145°F)
Shock	MIL-STD-202, Method 213 (50 g, half sine, 11 msec)
Vibration	MIL-STD-202, Method 204 (10 to 2000 Hz at 10 g)
Media	P2 port Wetted materials; polyester housing, epoxy adhesive, silicon, borosilicate glass, and silicon-to-glass bond*
	P1 port Dry gases only

* Liquid media containing some highly ionic solutions could potentially neutralize the chip-to-glass tube bond.


LEADWIRE OPTION

140PC sensors can be furnished with 12-inch, color coded, 24 gage leadwires instead of printed circuit board terminals. To order a leadwire version add a "W" to the end of a catalog listing shown in the order guide. For example, a 142PC15G with leadwires becomes 142PC15GW.

Pressure Sensors


Absolute, Differential, Gage, Vacuum Gage/Amplified

MOUNTING DIMENSIONS (For reference only)

Dimensions shown apply to Differential and Absolute versions. Gage units are identical, except the P1 port is absent.

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