Drehspiegel

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With the help of the rotary mirror (476 40), the lens of $f = \text{approx. } 5 \text{ m}$ (460 12), the front-silvered mirror (463 20), and some other accessories, the velocity of light can be measured according to the method devised by Foucault and Michelson.

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Notes

- The five-digit numbers quoted in parentheses refer to the catalogue numbers of the respective apparatus.
 - The data DC . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSICS LEAFLETS".
 Published to date:
- | | Cat. No. |
|----------------------------------|----------|
| 1st series (72 sheets) | 599 21 |
| 2nd series (52 sheets) | 599 22 |
| 3rd series (49 sheets) | 599 23 |
| 4th series (50 sheets) | 599 24 |

- The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the latest knowledge acquired in all scientific and technical fields.

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1. Description

The motor of the rotary mirror is housed in a casing of approx. 185 mm length and 55 mm \varnothing . The axle of the motor protrudes on top, and turns a mirror 10 mm \times 20 mm in size. The mirror is coated with an aluminium and quartz layer on either side. It is protected by a cap open on one side. The cap is provided with an adapter piece with a 12 mm hole for mounting the apparatus on stands.

A key can be used to adjust the mirror in the desired direction. The key is attached to the case with a chain and can be mounted on the mirror through a hole in the top of the cap.

The switch is at the bottom of the casing. The motor is permanently fitted with a two-core connecting cable that is approx. 7 m long and ends in Leybold plugs. The rotary mirror (476 40) can be connected to the 220 V mains through a 1000 Ω rheostat (537 22); rotary mirror (476 41) can be connected to 110 V mains through a 320 Ω rheostat (537 23). Both A.C. and D.C. are suitable for operation.

When the rheostat is not in circuit, the motor attains a maximum speed of approx. 28 000 r.p.m. It is provided with a lubricant that will last for the life of the motor and relubrication will not, therefore, be necessary.

The lens $f =$ approx. 5 m has a diameter of 120 mm and is fitted in a frame of approx. 180 mm external diameter. This frame is integral with a rod of 60 mm length and 10 mm \varnothing .

The front-silvered mirror has just the same frame as lens (460 12). It must be of particularly good optical quality in order to ensure good success of the experiment. The curvatures of the mirror surface, if any, are so slight as to result in focal lengths of not less than 100 m.

2. Arrangement

The experiment described in the following can be demonstrated in two alternative ways:

Without reflecting mirror, a straight light path of approx. 15 m length is required. An additional reflecting mirror reduces the linear extension of the arrangement to 7.5 m and also considerably facilitates the adjustment.

In addition to the apparatus mentioned above, the following pieces of equipment will be required:

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And additionally, for an arrangement
with reflecting mirror:

1 Stand base, 28 cm	300 01
1 Pair of levelling screws	300 06
1 Stand rod, 25 cm long	300 41
1 Leybold boss	301 01
1 Front-silvered mirror, $\varnothing 120$ mm, as reflecting mirror	463 20

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a) The projecting lens ($f = \text{approx. } 5 \text{ m}$) and the front-silvered mirror should each be attached to a 25 cm long stand rod which in turn should be inserted in the apex bore of the stand base. The stand rods, protruding somewhat out of the bore in the stand base, should serve as a point of support. The other two points of support are provided by the levelling screws for the stand base. The motor with rotary mirror should be fixed to a 50 cm long stand rod in the same manner.

b,) For the arrangement without reflecting mirror the apparatus should be arranged on three benches in the following manner:

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Bench 1:

Small optical bench in bench clamp with source of light, lens $f = 100$ mm, and adjustable slit. The lamp house should be pulled out as far as possible, and the lamp filament should be arranged in a vertical position and imaged on a wall nearly 5 m away by means of the lens. Directly behind the lens, provide a boss for attaching the slit later on. The lens $f =$ approx. 5 m should be set up to the side of this boss.

Bench 2:

Rotary mirror, approx. 5 m away from the boss for the slit on bench 1.

Bench 3:

Front-surface mirror as distant mirror, approx. 10 m away from the projecting lens.

b₂) For the arrangement with reflecting mirror the apparatus should be set up as follows:

Bench 1:
as under **b₁)**.

Bench 2:
Rotary mirror as in **b₁)**. The distant mirror should be set up beside the rotary mirror.

Bench 3:
Reflecting mirror, approx. 2.5 m away from the projecting lens.

c) The first stage in the adjustment of the optical system is to set the distance of the rotary mirror from the projecting lens by an autocollimation arrangement: The rotary mirror D should be illuminated by means of the lamp and condenser lens $f = 100$ mm and then turned with the key so that the reflected light falls on the centre of the projecting lens L set up about 5 metres away. The front-surface mirror should be placed a short distance behind this projecting lens so that the light beam reflected from the front-surface mirror forms an image of the rotary mirror alongside the rotary mirror itself. Find the location of the sharp image (holding a pin in front of the rotary mirror) and vary the distance between rotary mirror and lens L until image and rotary mirror are at equal distance from the lens. The correct positions for the rotary mirror and the lens will then have been fixed.

d₁) For the arrangement without reflecting mirror the above procedure should be followed by the adjustment of the distant mirror E₁ (see Fig. 2, full lines). At the selected distance $r = f$, the mirror must be two focal lengths f away from lens L. For checking the correct distance, the slit must briefly be re-inserted. When the

Table 1:

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light from the rotary mirror falls on the centre of the projecting lens, then a sharp image of the slit must appear on the centre of the distant mirror. (Check by holding white paper and pin in front of slit, displace distant mirror if necessary. It may be necessary to correct the position of the centering screws of the lamp house to obtain a better illumination). This having been done, manipulate the levelling screws of the stand base to bring the plane of the mirror in a position where the reflected light also passes through the centre of the lens. This adjustment can be facilitated considerably by holding a sheet of white paper horizontally in front of the lens.

d₂) With reflecting mirror: As shown in dashed lines in the assembly illustrated in Fig. 2, the reflecting mirror U should be arranged 2.5 m behind the lens, and the distant mirror E₂ should be set up alongside the rotary mirror D, in such a manner that the beam of light by-passes the projecting lens on its way to the distant mirror. The conditions that must be satisfied by the arrangement are the same as those for the arrangement without reflecting mirror: Light from the rotary mirror must fall on the centre of the projecting lens, then on the centre of the reflecting mirror U, and subsequently on the centre of the distant mirror E₂. To check the adjustment, re-insert the slit which must form a sharp image on the distant mirror. The definition (check with pin held before slit) can be improved by moving the reflecting mirror. It may be necessary to improve the illumination of the slit by shifting the lamp in the lamp house. The adjustment of the return path of the light where the image of the rotary mirror must be located on the rotary mirror is considerably easier to make because of the close proximity of rotary mirror and distant mirror.

e) The position of the distant mirror should now be checked on the rotary mirror, and finely adjusted. The image of the rotary mirror must be received on the rotary mirror itself. This image can be found by inserting a bit of paper in the beam of light directly in front of the rotary mirror. Slight lateral and vertical displacements of the distant mirror will bring the mirror and its image to exact coincidence.

f) The next stage is to set up the glass rule and the eyepiece. Care must be taken not to change the position of the stand bases and the various component parts.

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Insert the adjustable slit and arrange a clean glass plate (e. g. microscope slide) about 10 cm to 15 cm in front of it inclined at an angle of approx. 45° , using the holder with spring clips to fasten it on the small optical bench. To the rod of the holder with spring clips attach a 50 cm long stand rod and on this fix a second holder with spring clips into which the glass rule has been inserted. The glass rule should be in the beam of light reflected by the obliquely arranged glass plate in the position where the image of the slit is most well-defined. With a correct adjustment, therefore, it should be at the same distance from the oblique glass plate as the slit itself. The image of the now approx. 0.5 mm wide slit will be seen doubled. It should be examined through a second lens $f = 100$ mm used as an eyepiece. An iris diaphragm at the position of the aperture (image of the rotary mirror, approx. 100 mm behind the eyepiece) will facilitate bringing the eye into the right position for observation and also prevent the observation of undesirable scattered light. The glass rule should be so adjusted that the centre of the double image reflected at the object slide coincides with a graduation mark.

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g) The last stage is to connect the motor to the mains supply. Before doing so, make sure that the key on the mirror has been taken off. The motor should be connected to the A.C. mains in the manner shown in the circuit diagram (Fig. 3). On actuating switch S_1 , the motor will run slowly due to the series rheostat. At this speed, the image of the slit should be seen in the eyepiece at the same position as before, although it will be somewhat less bright. On throwing switch S_2 , the speed of the motor depends upon the position of the slider.

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h) The correct speed of the motor can be set with the help of an acoustic method involving the comparison of the noise of the motor with the sound of a 440 c/s tuning fork and adjusting the motor speed for zero beat. At this speed the deflection s should then be measured. It is also necessary to measure the focal length f , the lens not having a focal length of exactly 5 m.

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$$c = \frac{24f^2 n \pi}{s} \text{ m/s}$$

n = Drehzahl, U/s

The velocity of light can then be computed from the formula:

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3. Note

Be careful not to touch the front-surface mirror with the fingers, as this may result in the front coating being damaged. Any finger marks that might nevertheless be found on it should be removed very carefully with Green's lens tissue (30500).

The motor must never be run with the key inserted.

4. Experiment

Determination of the velocity of light according to Foucault and Michelson.

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