KIPPA d  
ZOVEN

Física

P18-28

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Técnica

LEYBOLD DIDACTIC GMBH

1/92 - rs -


Fig. 1

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Moll's thermopile is used to measure radiant heat energy. It consists of several series-connected thermocouples. When the pile is exposed to radiant heat energy, a thermoelectric voltage is built up which is measured by a voltage-sensitive measuring instrument.

**1 Technical Data**

Internal resistance: approx. 10 Ω  
 Sensitivity: approx. 0.16 mV/mW  
 Wavelength range: 150 nm – 15 µm  
 Response time: 2 to 3 secs.  
 Thermocouples: 16 thermocouples, series connected, composed of constantan-manganin strips, silver-soldered. Diameter of the measuring surface: 10 mm. Distance *a* (between thermocouples and front opening): approx. 50 mm (see Fig. 2)  
 Dimensions: Housing: 80 mm length, 34 mm dia.  
 Rod: 100 mm length, 10 mm dia.  
 Weight: approx. 0.5 kg

**2 Description**

**2.1 Functional parts (see Fig. 2)**

- ① Massive metal housing with rod
- ② Polished metal funnel with a 22° angle of aperture, min. transmittance aperture 10 mm
- ③ Blackened disk on the rear of which the soldered test junctions ⑤ of the 16 series-connected thermocouples ④ are positioned on the vertical plane diameter
- ④ Thermocouples


Fig. 2

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- ⑤ Soldered test junctions (thermocouple pairs soldered with silver)
- ⑥ Lacquer-coated plastic rods; used as heat conductors between thermocouple ④ and reference junction ⑦
- ⑦ Reference junction; is in direct thermal contact with housing ①
- ⑧ Electrical leads to measuring output
- ⑨ 4-mm sockets for connecting a sensitive measuring instrument
- ⑩ Glass disk in metal frame; is used as a protection during transport and storage and must be removed when performing measurements.

## 2.2 Operation (see Fig. 2)

The radiation incident on the thermopile strikes the blackened disk ③ partially direct and partially after deflection by the polished metal funnel ②. The blackened disk is in direct heat contact with the soldered test junctions ⑤ of the 16 series-connected thermocouples ④. The reference junction ⑦ is at the temperature of housing ①, i. e. virtually at room temperature.

If no radiant heat energy reaches the thermopile, soldered test junction ⑤ and reference junction ⑦ are at the same temperature. There is no temperature difference and hence no thermoelectric potential.

When now exposing the thermopile to the radiant heat energy of a body whose temperature is higher than that of reference junction ⑦ (housing or room temperature respectively), the blackened disk ③ rapidly becomes warm (in approx. 2 to 3 secs.) due to its low thermal capacity, while the housing remains at room temperature. This causes a temperature difference between soldered test junction ⑤ and reference junction ⑦. This results in a thermoelectric potential which can be tapped at the output sockets ⑨. The measured thermoelectric potential is proportional to the incident heat radiation power.

When now exposing the thermopile to the heat radiation of a body whose temperature is below the temperature of reference junction ⑦, the temperature of the blackened disk ③ will become lower than that of the reference junction (or housing respectively). Again there is a temperature difference from which now a thermoelectric potential of opposite polarity results.

## 3 Notes on Experiments

Suitable measuring instruments for the thermopile are:

- Demo - Multimeter (531 91)
- or
- Microvoltmeter (522 13)
- or
- Voltage sensitive measuring amplifier (522 06, no longer deliverable in conjunction with 60 mV/300µA movement, e.g. moving-coil instrument P (531 80)
- or
- School mirror galvanometer (532 10) in conjunction with the sensitivity regulator (532 11).

When using one of these instruments, please read the respective Instruction Sheets. Furthermore, observe the maximum load of the thermopile when performing experiments (see Technical Data). For precise measurements foreign influences (air draft or foreign heat radiation, such as emitted by sun or heater) should be avoided.