Planck'sche-Konstante-Apparat U10700

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05/10 ALF

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Instruction manual

05/10 ALF

- 1 LED with connecting lead
- 2 Nanoammeter
- 3 Voltmeter
- 4 Protective cover for photocell
- 5 Collector tube for photocell
- 6 Coaxial power socket for power supply
- 7 Power supply with connector socket for LED
- 8 Back EMF adjustment (coarse)
- 9 Back EMF adjustment (fine)
- 10 Intensity knob
- 11 Empty sleeve

1. Notes

The photocell in the Planck's constant apparatus is highly sensitive. Bright light can cause it to age quickly and can cause permanent damage. Following irradiation it ideally needs to be left for quite a long period before it regains its stability.

- The protective cover for the photocell should never be removed.
- When the experiment is completed, slide the empty sleeve over the collector tube of the photocell.
- Keep the apparatus secure so that it does not get shaken and do not expose it to extreme temperatures, high humidity, moisture or direct sunlight.

2. Scope of delivery

- 1 Basic apparatus with photocell, voltmeter, nanoammeter and LED power supply
- 1 Empty sleeve for covering the photocell collector tube
- 5 LEDs (472 nm, 505 nm, 525 nm, 588 nm, 611 nm) in case with connector leads
- 1 Plug-in power supply, 12 V AC
- 1 Instruction manual

3. Description

The Planck's constant apparatus is for determining the magnitude of Planck's constant and the work done in emitting electrons from a caesium cathode in a photocell using the back-EMF method.

It contains a vacuum photocell, a voltmeter for measuring back EMF, a nanoammeter for measuring the photocell current and a power supply for the LEDs. Five different light-emitting diodes (LEDs) are provided, which emit light at differing known average frequencies. The intensity of the emitted light can be varied between 0 and 100% in each case. The photocell itself consists of a cathode with caesium condensed onto its surface and a ring-shaped anode. When the apparatus is switched on, a voltage is applied between the two electrodes and this can be adjusted by two knobs for coarse and fine adjustment.

Power is supplied to the apparatus via the plug-in transformer provided. The Planck's constant apparatus with order number U10700-115 is designed for a mains voltage of 115 V ($\pm 10\%$) while the version with order number U10700-230 is for 230 V ($\pm 10\%$).

4. Technical data

Photocell:	Type 1P39, caesium (Cs)
Voltmeter:	3½-digit LCD
Precision:	0.5% (typically)
Nanoammeter:	3½-digit LCD
Precision:	1% (typically)
LEDs:	472 nm, 505 nm, 525 nm, 588 nm, 611 nm
Dimensions:	280x150x130 mm
Weight:	1.3 kg approx.

5. Theoretical principles

The light from the LED in the circuit passes through the ring-shaped anode and hits the cathode, where it causes emission of electrons with kinetic energy

$$E_{\text{kin}} = h \cdot f - W$$

According to Albert Einstein's photoelectric theory, $E = h \cdot f$ is the energy associated with photons of light at the wavelength in use.

W is the work done in emitting electrons from the cathode, i.e. the minimum energy required to expel electrons from the surface of the metal. This value is dependent on the material as well as on the temperature. For caesium it is 2.14 eV at 0 K and about 2 eV at room temperature.

Depending on the adjustment of the back EMF between the cathode and anode, a current of electrons should flow from the former to the latter. This can be measured using the nanoammeter. If the back-EMF corresponds to the critical voltage U_0 , where

$$e \cdot U_0 = E_{\text{kin}} = h \cdot f - W \quad \text{and } e = 1.6021 \times 10^{-19} \text{ C,}$$

then this current should have a magnitude of 0 nA.

Plot a graph of $e \cdot U_0$ against f for the critical voltages U_0 , measured for various frequencies of light f , to obtain a line of gradient h crossing the y axis at $-W$.

6. Operation

6.1 Measurement of critical voltage at a light intensity of 75%.

- Plug in the transformer to supply power.
- Set the intensity of the light source to 75%.
- Insert the plug for the first light source into the LED connector socket.
- Push together the jaws of the clip for the sleeve over the collector tube of the photocell and remove the sleeve.
- Push the LED unit fully onto the collector tube of the photocell until the jaws of the clip snap into place.
- Set the fine adjustment knob for the back-EMF to a central position.
- Wait a few minutes, then set the photocell current roughly to zero using the coarse adjustment knob.
- Find an optimum zero point using the fine adjustment knob.
- Take note of the back-EMF as set in this fashion and record it as the critical voltage U_0 .
- Repeat this measurement for the four other LEDs.

6.2 Determining Planck's constant h .

- Work out the frequencies of the light from the printed wavelengths λ using the formula $f = \frac{c}{\lambda}$

$$\text{where } c = 3 \cdot 10^8 \frac{\text{m}}{\text{s}}.$$

- Use the critical voltages U_0 to work out the energies $e \cdot U_0$ where $e = 1.6021 \times 10^{-19} \text{ C}$.
- Plot the values obtained on a graph of energy against frequency.
- Draw a straight line through the points and determine Planck's constant h from the gradient and the work W from where the line crosses the Y axis.

6.3 Proof that the critical voltage does not depend on the light intensity.

- Select an LED.
- Set the light to maximum intensity and determine the critical voltage U_0 .
- Reduce the intensity to zero in a set of steps and determine the critical voltage U_0 in each case.

7. Cleaning

- Only use conventional mild washing-up liquid and never aggressive detergents.
- Make very sure that no moisture gets inside the equipment.
- Unplug the power supply to turn off the equipment.
- Clean the equipment using a slightly damp and non-fluffy cloth.

8. Changing fuses

- Any blown fuses may only be replaced with a fuse matching the original.


Fig.1 Critical energy $U_0 \cdot e$ against frequency f


Fig. 2 Critical voltage U_0 against light intensity I at a wavelength of 472 nm

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