P43.5

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High tension power unit

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1. Description

The high tension power unit meets the requirements for a reliable source of variable high-tension D.C. for carrying out electrostatic experiments, for operating small gas discharge tubes, for post-acceleration in the cathode ray oscillograph for instructional use (575 15), and for other uses. It is operated from the A.C. mains, 220 volts, 50 cycles, and supplies two equal series-connected D.C. voltages which at the same time can be regulated continuously from 0 to approximately 3000 volts each with a rotary potentiometer. Built-in safety resistors limit the offtake of D.C. so that the high tension side can be short-circuited without causing any damage. A large band-pass filter renders the D.C. voltage practically free from humming.

The high tension power unit is located in a casing measuring 170 x 220 x 190 mm fitted at the top with a carrier handle. The front plate is fitted with the three high tension sockets, the rotary knob of the po-

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Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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potentiometer, the mains switch and an earth socket for the casing which is also earthed via the protective contact of the power supply plug. The high tension itself has no definite potential to earth. Any one of the three high tension sockets can be earthed.

2. Technical data

The basic circuit diagram of the high tension power unit is illustrated in Fig. 2.

Offtake outputs:

No-load voltage 0 ... 2 x approx. 3,000 volts
Short-circuit current approx. 2.5 milliamps.

Fuse:

Microfuse, rated value of 500 milliamps., 20 mm long, in rear plate of casing.

3. Note

The casing, which actually is free of voltage, must be earthed to protect whoever is carrying out the experiment. Then the earth socket of the front plate must be earthed reliably where it is in full view.

Owing to the large band-pass filter of the unit, the high tension sockets still carry high voltage even after the unit has been switched off. This high voltage does not fade away until several minutes have elapsed. If necessary the discharge should be accelerated by connecting the three high tension sockets in series. Owing to the installation of the large band-pass filter, the high voltage only decreases gradually when the knob is turned down on the rotary potentiometer.

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