# ORITEL

- PROPAGATION EN ESPACE LIBRE
- PROPAGATION GUIDEE
- PROPAGATION IN FREE SPACE
- WAVEGUIDE PROPAGATION

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# **H S - F G Z M**

Significance of symbol: 

**CAUTION! Refer to the experimentation manual before using the microwave instruction bench.**

In these operating instructions, failure to follow or carry out instructions preceded by this symbol may result in personal injury or damage to the device and the installations.

You have just acquired an **ORITEL microwave instruction system** and we thank you for your confidence.

To obtain the best possible service from this system :

- read this experimentation manual carefully,
- comply with the precautions for use.

## **PRECAUTIONS FOR USE**

- The ferrite isolator must be assembled correctly with the arrow pointing towards the load.
- Ferrite isolators and circulators are components with a magnetic core. Therefore, it is recommended to keep them at least 3 cm away from any magnetized object or source.
- When you connect the GUNN oscillator and the PIN Modulator to the power supply, any inversion of the power supply output connectors may damage the diodes (GUNN and PIN).
- The leads and connection accessories must comply with the applicable standards and must be designed for an overvoltage category rated voltage at least equal to that of the circuits on which the measurements are made.
- All the equipment connected to the bench and to the AC network must be earthed.
- No voltages of more than 30 V with respect to earth will be used on the bench.
- Although the transmission power is low, it is highly inadvisable to look directly, front on at operating transmission antennas

## **GUARANTEE**

Our guarantee is applicable, unless otherwise stated, for **twelve months** following the date of supply of the equipment (extract from our General Sales Conditions, available on request).

## SOMMAIRE

<b>1. ACTIVATION .....</b>	<b>66</b>
<b>1.1. BENCH SET-UP .....</b>	<b>66</b>
1.1.1 Use of EASYFIX <sup>(TM)</sup> quick attachment.....	66
<b>1.2. COMMISSIONING AND USE OF CF 204 POWER SUPPLY OR     MEASUREMENT EQUIPMENTS (OPTION) .....</b>	<b>66</b>
<b>2. GENERALITIES .....</b>	<b>67</b>
<b>2.1. MICROWAVE COMPONENTS.....</b>	<b>67</b>
2.1.1 ADZ 100/3 impedance matching device with 3 plungers.....	67
2.1.2 LAZ 100 slide impedance adapter.....	67
2.1.3 ATC 100 calibrated attenuator.....	67
2.1.4 ATM 100 micrometer variable attenuator .....	67
2.1.5 CAB 100 coaxial cable.....	67
2.1.6 CHG 100 matched load.....	67
2.1.7 COE 100/H/B E plane bend (high and low .....	67
2.1.8 COH 100 H plane bend .....	67
2.1.9 CDT 100/20 configurable hole directive coupler .....	67
2.1.10 CGX 100/20 configurable cross directive coupler.....	67
2.1.11 IRIS 20.30 cross coupling 20 and 30 dB iris .....	67
2.1.12 ICDT 100/30 30 dB hole coupling irises.....	67
2.1.13 CIR 100 ferrite circulator .....	67
2.1.14 ISO 100 ferrite isolator.....	68
2.1.15 CC 100 short-circuit.....	68
2.1.16 CCM 100 micrometer variable short-circuit .....	68
2.1.17 DPH 100 micrometer phase-shifter .....	68
2.1.18 DEN 100 coaxial detector.....	68
2.1.19 DEG 100 waveguide parallel detector.....	68
2.1.20 AFR 100 fastening adapter .....	68
2.1.21 GD 100/180 straight waveguide.....	68
2.1.22 JTG 100 rotary joint .....	68
2.1.23 KED 100 dielectric study kit.....	68
2.1.24 LAF 100 slotted line .....	68
2.1.25 MOD 100 PIN diode modulator .....	68
2.1.26 OND 100 micrometer wavemeter.....	69
2.1.27 OSG 100 Gunn diode oscillator.....	69
2.1.28 TGN 100 waveguide-coaxial transition .....	69
<b>2.2. MICROWAVE ANTENNAE .....</b>	<b>69</b>
2.2.1 AND 100 dielectric antenna .....	69
2.2.2 ANF 100 slot antenna .....	69
2.2.3 ANF 100F slot antenna .....	69
2.2.4 ASP 100 microstrip antenna .....	69
2.2.5 ANC 100/10 (or ANC 100/15 or ANC 100/20) horn antennae.....	69
2.2.6 DR 100 reflector disk .....	69
2.2.7 ANP 100 parabolic reflector.....	69
2.2.8 ANP 100 F - Fixed parabolic reflector .....	69
2.2.9 RRL 100 passive RADAR responder.....	69
2.2.10 Adjustable slot iris.....	69
<b>2.3. SUPPORTS AND ACCESSORIES .....</b>	<b>70</b>
2.3.1 BTE 100 experimentation frame .....	70
2.3.2 ABS 100 absorbing panels .....	70
2.3.3 PTM 100 manual rotating platform.....	70
2.3.4 SAN 100 antenna support .....	70
2.3.5 SUP 100 waveguide support.....	70
<b>2.4. POWER SUPPLY AND MEASURING INSTRUMENTS.....</b>	<b>70</b>

2.4.1	CF 204 - GUNN power supply .....	70
2.4.2	IR 205 SWR indicator.....	70
2.4.3	MH 600 microwave powermeter .....	70
2.4.4	ST 613 thermocouple coaxial probe (to be used with the MH 600 powermeter) .....	70
2.4.5	SDD 100 displacement simulator.....	71
2.4.6	Oscilloscopes: Analogue and digital: OX 520BS – MTX 3252M – MTX 3352C .....	71
<b>3.</b>	<b>GENERAL .....</b>	<b>72</b>
3.1.	BENCH INSTALLATION.....	72
3.2.	BEM 100 TRANSMISSION BENCH .....	73
<b>4.</b>	<b>EXPERIMENTATION.....</b>	<b>74</b>
4.1.	SPACE ATTENUATION MEASUREMENT.....	74
4.1.1	Introduction.....	74
4.1.2	Theoretical reminders .....	74
4.1.3	Space attenuation study .....	74
4.2.	POLARIZATION OF AN ELECTROMAGNETIC WAVE.....	75
4.2.1	Introduction.....	75
4.2.2	Theoretical reminders .....	75
4.2.3	Polarization study.....	75
4.3.	GAIN MEASUREMENTS .....	77
4.3.1	Introduction.....	77
4.3.2	Theoretical reminders .....	77
4.3.3	ANTENNA GAIN MEASUREMENT (Method using two identical antennae).....	79
4.4.	HORN ANTENNA STUDY .....	81
4.4.1	Introduction.....	81
4.4.2	Theoretical reminders .....	81
4.4.3	Radiation pattern measurement.....	81
4.4.4	Gain Measurement.....	82
4.5.	DIELECTRIC ANTENNA STUDY .....	83
4.5.1	Introduction.....	83
4.5.2	Theoretical reminders .....	83
4.5.3	Radiation pattern measurement.....	84
4.5.4	Gain measurement .....	85
4.6.	STUDY OF SLOT ANTENNA .....	86
4.6.1	Introduction.....	86
4.6.2	Theoretical reminders .....	86
4.6.3	Babinet principle.....	89
4.6.4	Radiation pattern measurement.....	90
4.6.5	Gain measurement .....	93
4.7.	PARABOLIC ANTENNA STUDY.....	94
4.7.1	Introduction.....	94
4.7.2	Theoretical reminders .....	94
4.7.3	Radiation pattern measurement.....	96
4.7.4	Gain measurement .....	98
4.8.	PASSIVE RADAR RESPONDER STUDY .....	99
4.8.1	Introduction.....	99
4.8.2	Theoretical reminders .....	99
4.8.3	Measurement of RSE (Equivalent Radar Surface) .....	102
4.9.	STUDY OF DOPPLER EFFECT .....	104
4.9.1	Introduction.....	104
4.9.2	Theoretical reminders .....	104
4.9.3	Speed of displacement measurement.....	106
4.10.	DIRECTIVE COUPLER STUDY.....	108
4.10.1	Introduction.....	108
4.10.2	Theoretical reminders .....	108
4.10.3	Coupler experimentation .....	110
4.11.	STUDY OF FERRITE ISOLATOR AND CIRCULATOR .....	114
4.11.1	Introduction.....	114
4.11.2	Theoretical reminders .....	114
4.11.3	Experimentation.....	115

# INTRODUCTION

Because of an increasing demand for microwave technique training, the ORITEL division at CHAUVIN ARNOUX has developed measurement equipment specially designed for training purposes.

ORITEL, having specialized in microwave components or systems for several years, was concerned to simplify the various components used in these measurement benches. This set of techniques is therefore now within reach of a wide range of users.

The frequency band: 8.5 GHz to 9.6 GHz ( $\lambda = 3$  cm) was chosen first because of the major developments in this frequency range and then, because of the size of the components.

Other components available from the CHAUVIN ARNOUX catalogue can be supplied on demand.

Far from being a microwave measurement course, the sole purpose of this experimentation manual is to illustrate, through a few examples, the possibilities of the ORITEL equipment.

Note that this training bench has one innovating particularity which greatly facilitates the assembly of its components: the **EASYFIX™** fast attachment system which, by a single manoeuvre, makes it possible to connect any microwave components to form a waveguide by coupling connecting flanges of the UBR100 / UG39 (R100/WR90) type together.

Fast attachment  
adapter EASYFIX™

# 1. ACTIVATION


## 1.1. BENCH SET-UP

Because the assembly of the various components differs, depending on the measurements to be made, details will be given for each of the manipulations described in this manual.

As a general rule, whatever the measurement, it is mandatory to install the ferrite isolator immediately after the GUNN oscillator so that the frequency and power stabilities of the oscillator are not or are only slightly affected by any impedance placed after the isolator.

### 1.1.1 Use of EASYFIX<sup>(TM)</sup> quick attachment

- Arrange the two parts to be assembled, by aligning correctly centring studs ① in front of the holes ②
- Simultaneously tighten the two thumbwheels ③

## 1.2. COMMISSIONING AND USE OF CF 204 POWER SUPPLY OR MEASUREMENT EQUIPMENTS (OPTION)

For the commissioning of the GUNN CF 204 power supply, the IR 205 VSWR indicator, the MH 600 powermeter and the OX 0520BS, MTX 3252M, and MTX 3352C oscilloscopes, refers to the corresponding operation manuals.

## 2. GENERALITIES

### 2.1. MICROWAVE COMPONENTS

#### 2.1.1 ADZ 100/3 impedance matching device with 3 plungers

This matching device consists of a waveguide section with 3 screw plungers to generate impedance mismatch in the frequency band.

#### 2.1.2 LAZ 100 slide impedance adapter

Inserted into a microwave line, this component matches the impedance of any load and, accordingly, decreases the standing wave ratio.

#### 2.1.3 ATC 100 calibrated attenuator

This variable calibrated attenuator is designed to proportion the energy supplied by the oscillator. The attenuation induced is indicated by means of a dial graduated in dB.

#### 2.1.4 ATM 100 micrometer variable attenuator

The variable attenuator mounted on the waveguide is used for proportioning the energy transmitted over the line and making comparative level measurements in dB. The micrometer adjustment level variation control provides high accuracy and good adjustment repetition.

#### 2.1.5 CAB 100 coaxial cable

This complementary cable is used for guided propagation measurements in a coaxial line.

#### 2.1.6 CHG 100 matched load

This matched load consisting of a closed straight waveguide section containing a matched dissipating element, terminates the line on its characteristic impedance.

#### 2.1.7 COE 100/H/B E plane bend (high and low

Used for a changing the direction of the microwave path in the E plane.

#### 2.1.8 COH 100 H plane bend

Used for a changing the direction of the microwave path in the H plane.

#### 2.1.9 CDT 100/20 configurable hole directive coupler

This coupler has the particularity of having interchangeable coupling irises (20 and 30 dB).

#### 2.1.10 CGX 100/20 configurable cross directive coupler

This coupler has the particularity of having interchangeable coupling irises (20 and 30 dB).

#### 2.1.11 IRIS 20.30 cross coupling 20 and 30 dB iris

Coupling iris for CGX 100 coupler.

#### 2.1.12 ICDT 100/30 30 dB hole coupling irises

Coupling iris for CDT 100 coupler

#### 2.1.13 CIR 100 ferrite circulator

This element is a hybrid junction with the property of switching microwave energy in a privileged direction.

#### **2.1.14 ISO 100 ferrite isolator**

The ferrite isolator protects the GUNN diode oscillator from load impedance variations and from the mismatch caused by the PIN diode modulator.

The ferrite isolator allows microwave energy to pass without attenuation in the oscillator-load direction and generates strong attenuation in the opposite direction.

#### **2.1.15 CC 100 short-circuit**

This board generates a perfect short circuit in the plane of reference of the transmission line consisting of waveguide R 100.

#### **2.1.16 CCM 100 micrometer variable short-circuit**

This element is intended to move the reference plan of the short circuit along the line of transmission, in order to generate adaptations of impedance, phase shift and reflected waves.

#### **2.1.17 DPH 100 micrometer phase-shifter**

This component is designed to introduce a phase-shift to a transmission line. A micrometer is used for adjusting the value.

#### **2.1.18 DEN 100 coaxial detector**

This detector is used to detect the microwave energy. A Schottky diode ensures amplitude detection while the BNC coaxial plug it is used for connecting the device measuring the detected voltage. Mounted on the ORITEL TGN 100 waveguide-coaxial transition, it establishes a detector on the waveguide R100.

#### **2.1.19 DEG 100 waveguide parallel detector**

This detector is designed for studies of the Doppler effect. It is used for revealing the cross-modulation products required for this study.

#### **2.1.20 AFR 100 fastening adapter**

Equipped with the quick attachment system EASYFIX<sup>(TM)</sup> this module allows the quick assembly of 2 components consisting of R100 / WR90 waveguide sections terminating in the connection flanges UBR100/UG39.

#### **2.1.21 GD 100/180 straight waveguide**

This straight waveguide section, 180 mm in length, is indispensable for the study of dielectric materials.

#### **2.1.22 JTG 100 rotary joint**

Necessary for measuring an antenna radiation pattern, it allows antenna rotation in the horizontal plane through  $\pm 180^\circ$  (or  $360^\circ$ ).

#### **2.1.23 KED 100 dielectric study kit**

The purpose of the accessories is to study the effect of the dielectric environment on the guided propagation properties. It consists of samples of different dielectric materials of a variety of shapes.

#### **2.1.24 LAF 100 slotted line**

The slotted line is used for measuring the standing wave ratio amplitude and making impedance measurements.

A detector probe inserted into the waveguide through a longitudinal slot picks up part of the energy and detects it. The carriage supporting allows the move over the entire length of the slot.

#### **2.1.25 MOD 100 PIN diode modulator**

The PIN diode modulator is designed to ensure amplitude modulation by square wave signals at 1kHz, of the microwave produced by the GUNN oscillator.

### **2.1.26 OND 100 micrometer wavemeter**

The absorption wavemeter consisting of a resonating cavity coupled with a waveguide by an iris ensures frequency measurement of the energy derived from the oscillator setup. Volume adjustment in the cavity is by a micrometer graduated in mm, while a calibration curve specifies the typical tuning frequency of the cavity according to the graduation.

### **2.1.27 OSG 100 Gunn diode oscillator**

The Gunn diode oscillator is a solid-state microwave energy source. This oscillator consists of a waveguide section one end of which is closed by a mobile short-circuit, being frequency-adjustable from 8.5 to 9.6 GHz by means of a micrometer control.

### **2.1.28 TGN 100 waveguide-coaxial transition**

The waveguide-coaxial transition is a way of terminating the waveguide with a N female type coaxial connector.

## **2.2. MICROWAVE ANTENNAE**

### **2.2.1 AND 100 dielectric antenna**

This antenna is a radiating waveguide made of dielectric material. By combining 2, 3 or 4 antennae together, AND 100 offers the possibility of constructing a dielectric antenna array.

### **2.2.2 ANF 100 slot antenna**

The radiating slots in this antenna can be used or masked to study variations in the radiation pattern and the gain, depending on the number of slots being used.

### **2.2.3 ANF 100F slot antenna**

This antenna is identical to ANF 100 but will not allow the slots to be masked.

### **2.2.4 ASP 100 microstrip antenna**

This antenna uses the "microstrip" technology, by the association of individual dipoles (patches).

### **2.2.5 ANC 100/10 (or ANC 100/15 or ANC 100/20) horn antennae**

These pyramidal antennae have respective gain values of 10,15 and 20 dB. These horns are used in a variety of setups for radiation pattern measurement on the transmission bench, propagation measurement in free space, parabolic antenna illumination, etc

### **2.2.6 DR 100 reflector disk**

This 12.7 mm diameter disk is used for experimentation on the passive radar responder, and for other manipulations requiring a reflector.

### **2.2.7 ANP 100 parabolic reflector**

This 350 mm diameter reflector is adjustable through approximately 15° in elevation. The position of the source around the focus point is also adjustable. Illumination is of the axial type.

### **2.2.8 ANP 100 F - Fixed parabolic reflector**

This reflector is identical to ANP 100 but will not allow elevation or focal length adjustments.

### **2.2.9 RRL 100 passive RADAR responder**

This responder has the property of reflecting the energy it receives with constant amplitude over a large angular sector. It is as if it concentrates the energy it receives, reflecting it in the direction of the incident wave, through an aperture angle of around 180°.

### **2.2.10 Adjustable slot iris**

Used for experimentation of ANF 100 using an ANF 100F antenna.

## 2.3. SUPPORTS AND ACCESSORIES

### 2.3.1 BTE 100 experimentation frame

This frame holds the antennae at a distance from surfaces that can cause interference by reflection. It includes shelves designed to support the measurement accessories.

### 2.3.2 ABS 100 absorbing panels

During antenna radiation measurements, these panels are used for protecting the receiver from unwanted reflection whose amplitude may equal that of the effective signal.

### 2.3.3 PTM 100 manual rotating platform

This manual rotating platform is designed to turn all the proposed antennae through  $\pm 180^\circ$  in order to study them; it is provided with an indispensable SAN 100 antenna support

### 2.3.4 SAN 100 antenna support

This device maintains the various antennae proposed to allow their study.

### 2.3.5 SUP 100 waveguide support

This device is designed to support the components of a bench when it is being used on a working surface.

## 2.4. POWER SUPPLY AND MEASURING INSTRUMENTS

### 2.4.1 CF 204 - GUNN power supply

This power supply is initially designed to power the Gunn effect diode oscillator OSG 100, and the P.I.N. diode modulator MOD 100.

It is specifically intended to avoid destruction of the Gunn diodes when they are powered up, or when characteristics are measured.

### 2.4.2 IR 205 SWR indicator

The SWR indicator is a selective voltmeter whose frequency is adjustable around 1000 Hz. This indicator can be used for four main functions:

Measurement of the Standing Wave Ratio (SWR) using a slotted line of the LAF 100 type

- Insertion loss measurement
- Measurement of an attenuation
- Indication of the field variations received by an antenna

### 2.4.3 MH 600 microwave powermeter

This microwave digital powermeter:

- Measures power over a dynamic range of 50 dB, in dBm or mW, and power deviations in dB or %, in the LF, HF, UHF and SHF bands,
- Uses the ORITEL ST 600 series thermocouple probes, the ORITEL SD 600 series Schottky diode probes.
- Is compatible with probes from the HP8480 family.

### 2.4.4 ST 613 thermocouple coaxial probe (to be used with the MH 600 powermeter)

The ST 613 probe is one of the microwave power measuring elements of the ORITEL MH 600 digital powermeter. This probe, measuring from 10  $\mu$ W to 100mW, between 10 MHz and 18 GHz, uses a microwave thermocouple using thin film technology.

#### **2.4.5 SDD 100 displacement simulator**

To study the Doppler effect, it is essential to have a moving target, whose speed of movement can be varied, as well as a large enough room to conduct of the experiment.

The movement simulator avoids the drawbacks related to the movement of an object over considerable distances and with variable speeds; this simulator is static in appearance.

It consists of a case containing a dielectric wheel provided with metal contacts, driven in rotation by a motor whose speed can be varied by a simple electric control.

The wheel is placed horizontally and parallel to the direction of microwave energy propagation; it is masked over approximately 80% of its height by an electrostatic screen so that only the upper part is in the electromagnetic field; in this way, the horizontal tangential rate of movement of the unmasked metal contacts simulates a linear speed of movement.

#### **2.4.6 Oscilloscopes: Analogue and digital: OX 520BS – MTX 3252M – MTX 3352C**

One of these oscilloscopes is indispensable for of the study of the Doppler effect.

### 3. GENERAL

#### 3.1. BENCH INSTALLATION

Essentially, installation consists in setting up the "Transmission bench + Reception bench" system in such a way as to establish a direct radiation link between them.


To minimize the effects of indirect radiation caused by reflection, the best arrangement consists in placing the transmission bench in an area clear of any obstacles which could form a source of reflection and enclosing it by absorbent panels on either side.

The distance R separating the transmission antenna from the reception antenna must be greater than  $2d^2/\lambda_0$

d = the largest dimension of the radiating aperture in the larger of the antennas.

$\lambda_0$  = the wavelength in the air.

A shorter distance R will induce an error in the measurement of the radiation patterns and in the gain measurements. Place the two assemblies so that they are correctly aligned with the antennas on exactly the same axis of symmetry. Level and stabilize the two experiment supports using the devices provided for the purpose.

The assembly is then ready to operate and to be used.

Proceed as follows:

- Ensure that the turntable rotates freely.
- Turn the reception system to obtain a maximum detected signal level on the SWR indicator (IR 205).

Adjust the "Frequency" setting of the indicator to obtain maximum of reflection on the galvanometer. Change the sensitivity if necessary.

If the detected signal becomes too strong, decrease the transmission level using an ATM 100 attenuator or the "level" control of the CF 204 power supply, for instance.

- Adjust the index of the turntable to zero so that it corresponds to the maximum deflection on the galvanometer.
- Adjust the gain setting of the SWR indicator to bring the galvanometer pointer to the end of the scale on the 0 dB graduation.
- When that these conditions are fulfilled, the measurement bench is ready to make measurements on the antenna.

**NOTE:** If the level differences are small, use the "extended dB" sensitivity scale on IR 205.

The measurements proposed in this manual are for the greater part presented using a coaxial detector (DEN 100) connected to an SWR indicator (IR 205). Note that other means of measurement can be used indeed. A powermeter (MH 600) equipped with a probe (ST 613) can be an advantageous replacement for the indicator (IR 205) and the coaxial detector (DEN 100). It is then necessary to transpose the measurement methods to this type of apparatus.

## 3.2. BEM 100 TRANSMISSION BENCH

### ■ Composition:

- 1 20 dB horn antenna: ANC 100/20
- 1 quick fastening adapter: AFR 100
- 1 PIN modulator: MOD 100
- 1 insulator: ISO 100
- 1 Gunn source: OSG 100
- 2 guide supports: SUP 100
- 1 GUNN power supply: CF204 (option)

### ■ Starting up the transmission bench

The transmission bench consists of an ANC 100/20 horn with an aperture of 163 x 114, pyramidal, powered by an OSG 100 Gunn diode oscillator. The signal is modulated by 1 kHz square signals by a MOD 100 pin diode modulator.

An ISO 100 ferrite insulator between the oscillator and the PIN via modulator prevents the load variations caused by the latter from reacting on the oscillator.

- Assemble the microwave components forming the BEM 100 transmission bench (fig. 1).
- Connect the OSG 100 Gunn diode oscillator from the CF 204 power supply by means of a BNCm/BNCm cord.
- Connect the MOD 100 Gunn diode oscillator from the CF 204 power supply by means of a BNCm/BNCm cord.
- Switch on the CF 204 power supply and increase the voltage gradually to 10V.
- The transmission bench transmits a square signal modulated microwave at a frequency of approximately 1 kHz by means of the ANC 100/20 horn towards the reception bench.


fig : 1

**Note:** To avoid unwanted reflection due to the nearness of the ground, which may interfere with the measurement, it is advisable to place the transmission and reception benches on BTE 100 experiment frames (option) and use absorbing panels to attenuator the spurious reflections.

## 4. EXPERIMENTATION

### FREE SPACE PROPAGATION

#### 4.1. SPACE ATTENUATION MEASUREMENT

##### 4.1.1 Introduction

The experiment consists in highlighting the attenuation brought in between transmission and reception according to the distance between two antennae in free space.

##### 4.1.2 Theoretical reminders

Referring to the formula:

$$Pr = \frac{Pe G^2 \lambda_0^2}{(4\pi)^2 R^2}$$

$\lambda_0$   $\lambda$  = wavelength in air.

G = Gain of each horn (identical)

We obtain:

$$\frac{Pr}{Pe} = \frac{G^2 \lambda_0^2}{(4\pi)^2} \cdot \frac{1}{R^2} = k \frac{1}{R^2}$$

K = Constant of experiment.

By varying R (distance between the two antennae), the  $Pr/Pe$  ratio will vary

The experiment consists in determining:

$$\frac{Pr}{Pe} \text{ for } R \quad ; \quad \frac{Pr_1}{Pe_1} \text{ for } R_1 \quad ; \quad \frac{Pr_2}{Pe_2} \text{ for } R_2 \quad ; \quad \text{etc...}$$

Transforming these results into dB and comparing the results found with the reading made on the calibrated attenuator (ATC 100).

##### 4.1.3 Space attenuation study

###### ■ Composition:

- 1 GUNN Source :	OSG 100	
- 1 Ferrite isolator :	ISO 100	
- 1 PIN Modulator :	MOD 100	
- 1 Micrometer attenuator :	ATM 100	
- 1 Micrometer wavemeter :	OND 100	
- 1 Calibrated attenuator :	ATC 100	
- 2 15 or 20 dB horn antennae :	ANC 100/15 or	ANC 100/20
- 1 Waveguide/coaxial transition:	TGN 100	
- 1 Capteur de mesure :	DEN 100	or 1 Thermocouple probe: ST 613
- 4 Support guide :	SUP 100	
- 2 Fastening adapters :	AFR 100	
- 1 Measurement indicator:	IR 205	or 1 Powermeter: MH 600
- 1 GUNN power supply:	CF 204	

■ **Activation of the bench**

The components used for the experiments that follow are essentially BDH R100 basic bench components.

- Assemble the microwave components and connect the CF 204 GUNN power supply and the IR 205 (or MH 600 and ST 613) measurement indicator as indicated below (fig. 1).
- Power up the GUNN power supply and the measurement indicator.
- Make the measurements.


fig 1

Compare the results obtained with the ATC 100 calibrated attenuator.

## 4.2. POLARIZATION OF AN ELECTROMAGNETIC WAVE

### 4.2.1 Introduction

The measurement of the attenuation due to the variation of electromagnetic wave polarization consists in determining the captured power variation according to rotation angle of the horn about its longitudinal axis of symmetry.

### 4.2.2 Theoretical reminders

The radiation polarization can be revealed by a simple application of the MALUS Law, well known to opticians.

The transmitted wave is polarized and rectilinear because the propagation mode in the rectangular waveguide is TE<sub>10</sub>. The received wave is rectilinear. If the receiver horn turns about its longitudinal axis through an angle  $\alpha$ , the power received will have the expression :

$$P = P_0 \cos^2 \alpha$$

$P_0$  being the power received without rotation of the receiver horn ( $\alpha = 0^\circ$ ).  
Attenuation then becomes:

$$A(\text{dB}) = 10 \log \frac{P_0}{P_0 \cos^2 \alpha} = -20 \log |\cos \alpha|$$

### 4.2.3 Polarization study

■ **Composition :**

- |                                    |                          |                                 |
|------------------------------------|--------------------------|---------------------------------|
| - 1 GUNN Source :                  | OSG 100                  |                                 |
| - 1 Ferrite isolator :             | ISO 100                  |                                 |
| - 1 PIN Modulator :                | MOD 100                  |                                 |
| - 1 Micrometer attenuator :        | ATM 100                  |                                 |
| - 1 Calibrated attenuator :        | ATC 100                  |                                 |
| - 2 15 or 20 dB horn antennae :    | ANC 100/15 or ANC 100/20 |                                 |
| - 1 E/H positionneur :             | PEH 100                  |                                 |
| - 1 Waveguide/coaxial transition : | TGN 100                  |                                 |
| - 1 Measurement sensor :           | DEN 100                  | or 1 Thermocouple probe: ST 613 |
| - 4 Waveguide supports :           | SUP 100                  |                                 |
| - 2 Fast attachment adapters :     | AFR 100                  |                                 |
| - 1 Measurement indicator :        | IR 205                   | or 1 Powermeter: MH 600         |
| - 1 GUNN power supply :            | CF 204                   |                                 |

■ **Activation of the bench**

The components used for the experiments that follow are essentially BDH R100 basic bench components.

- Assemble the microwave components and connect the CF204 GUNN power supply and the IR205 (or MH600 and ST613) measurement indicator as indicated below (fig. 1).
- Power up the GUNN power supply and the measurement indicator.
- Make the measurements.



fig : 1

Proceed as indicated in the figure below :

- Set the PEH 100 positioner to zero, corresponding stringently to the vertical axis of the horn.
- Set the ATC 100 attenuator to its zero position.
- Turn the receiver horn through an angle  $\alpha$  and calculate the corresponding attenuation.
- Locate the indication given on the IR 205 galvanometer.
- Bring the receiver horn to the position  $\alpha = 0$ .
- Adjust the ATC 100 attenuator to obtain the indication identified previously and locate the attenuator indication.
- Perform several manipulations, giving  $\alpha$  different values each time.

## 4.3. GAIN MEASUREMENTS

### 4.3.1 Introduction

The following measurements are designed to determine the gain value of an antenna using two different methods.

- Use of two identical antennae - as well as the attenuation of space.
- Three-antenna method.

### 4.3.2 Theoretical reminders

The antennae or aeriels are radiating devices ensuring propagation transition between the waveguide and free space (transmission) and vice versa (reception).

In the microwave field, because the wavelength is small, a great variety of radiating devices is possible: horns, parabolas, plane and dielectric antennae.

In reality, they are similar to optical devices: lenses, projectors for which, as in the microwave field, it is desired to concentrate energy in a given direction.

#### ■ Definition of antenna gain

If  $P_o$  is the power transmitted by a purely theoretical antenna which radiates uniformly in all directions (isotropic antenna), the gain of a radiated power antenna  $P(\theta, \phi)$  ( $\theta$  and  $\phi$  designating any direction in spherical coordinates) can be expressed :

$$G_{(\theta, \phi)} = \frac{P_{(\theta, \phi)}}{P_o} \quad (1)$$

The gain is a way of expressing the concentration of the energy transmitted (or received) in a privileged direction.

Therefore, an antenna is associated with its equivalent surface area (or effective surface area).

$$\Sigma(\theta, \phi) = G_{(\theta, \phi)} \frac{\lambda^2}{4 \pi} \quad (2)$$

For a theoretical isotropic antenna  $G_{(\theta, \phi)} = 1$  by definition, and therefore:

$$\Sigma(\theta, \phi) = \frac{\lambda^2}{4 \pi} \quad (3)$$

#### ■ Particular case of the horn antenna

The horn antenna shown in the figure opposite is obtained by the gradual dilatation of a conventional waveguide section, thus increasing the equivalent surface area and therefore the radiating aperture gain.

For a rectangle wave-guide of dimension  $a$  and  $b$  :

$$a < \lambda \text{ et } b \leq a / 2$$


The  $ab$  surface area is equal t :

$$ab \leq \lambda^2 / 2 .$$

Given this formula (2), the gain cannot exceed  $2 \pi$

Axial length  $l$  must also be such that the electric and magnetic field phase is constant in the horn output plane; this condition is obtained for:

$$l \geq d^2 / 2\lambda_0$$

#### ■ Measurement condition

In practice, two identical horns in terms of sizes are used, to reproduce remote field propagation conditions with spacing at a distance  $R$  such that :

$$R \geq 2d^2 / \lambda_0$$

$d$  : represents the largest size of the radiating apertures.

Although they are relatively bulky, to have high gain, the horns are easily adapted to the excitation waveguide and have a wider frequency band.

■ **Method using a pair of identical antennae**

If  $P_e$  is the transmitted power and  $P_r$  the received power,  $G$  the gain of the two identical antennae and  $R$  the separation distance, power  $P_r$  will be as given in the following relation :

$$P_r = \frac{P_e G^2 \lambda_0^2}{(4\pi)^2 R^2} \quad (4)$$

Immediately it is evident that :

$$G = \frac{4\pi R}{\lambda_0} \sqrt{P_r / P_e}$$

Or:

$$G_{(dB)} = \frac{1}{2} (P_{r(dB)} - P_{e(dB)}) + 10 \log (4\pi R / \lambda_0) \quad (5)$$

Furthermore, the expression (4) clearly reveals the dependence of  $1/R^2$  with respect to  $(P_r/P_e)$

$$\frac{P_r}{P_e} = (G \lambda_0 / 4\pi)^2 \frac{1}{R^2} \quad (6)$$

■ **Three-antennae Method**

The gain measurement methods generally proposed involve two identical antennae for transmission and reception and therefore mean that a pair of antennae has to be available. It is possible to, however, to make gain measurements on different antennae as long as there are three of them (we will refer to A, B and C).

Having made an initial measurement with antennae A and B, B is exchanged with C for a second measurement and finally, A with B for a third, while taking steps to maintain a constant wavelength and distance  $R$  between the antennae. This produces three equations with three unknowns :

$$G_A G_B = \left[ \frac{4\pi R}{\lambda_0} \right]^2 \left[ \frac{P_r}{P_e} \right]_1 \quad (1)$$

$$G_A G_C = \left[ \frac{4\pi R}{\lambda_0} \right]^2 \left[ \frac{P_r}{P_e} \right]_2 \quad (2)$$

$$G_C G_B = \left[ \frac{4\pi R}{\lambda_0} \right]^2 \left[ \frac{P_r}{P_e} \right]_3 \quad (3)$$

$\left[ \frac{P_r}{P_e} \right]_i$  being the ratio of the received and transmitted power in case  $i$ , with  $i = 1, 2$  and  $3$ .

the resolution of this system gives :

$$G_A = \left[ \frac{4\pi R}{\lambda_0} \right] \sqrt{\left[ \frac{P_r}{P_e} \right]_1 \left[ \frac{P_r}{P_e} \right]_2 \left[ \frac{P_e}{P_r} \right]_3}$$

i.e. in dB :

$$G_{A(dB)} = \left[ \frac{4\pi R}{\lambda_0} \right]_{(dB)} + \frac{1}{2} \left( \left[ \frac{P_r}{P_e} \right]_{1(dB)} + \left[ \frac{P_r}{P_e} \right]_{2(dB)} - \left[ \frac{P_e}{P_r} \right]_{3(dB)} \right)$$

We obtain  $G_B$  and  $G_C$  easily in the same way.

■ **Procedure using three different antennae (horn, slot and parabola):**

Make a gain measurement using, for instance, horn antenna (A) for transmission and the slot antenna (B) for reception. The antennae are placed 9 m apart, and it is essential to maintain this distance throughout the experiment.

Align the system by obtaining maximum deviation on the IR 205 level indicator.

Determine the ratio  $\left[ \frac{Pr}{Pe} \right]_1$  (A → B)

Disassemble the reception slot antenna and replace it by the parabolic antenna (C).

Work out the ratio  $\left[ \frac{Pr}{Pe} \right]_2$  (A → C)

Disassemble the transmission horn and replace it by the slot antenna.

Determine the ratio  $\left[ \frac{Pr}{Pe} \right]_3$  (B → C)

Calculate  $G_A$  with the formula given above.

Work out:

$$G_{B(dB)} = \left[ \frac{4\pi R}{\lambda_0} \right]_{(dB)} + \left[ \frac{Pr}{Pe} \right]_{1(dB)} - G_{A(dB)}$$

as well as :

$$G_{C(dB)} = \left[ \frac{4\pi R}{\lambda_0} \right]_{(dB)} + \left[ \frac{Pr}{Pe} \right]_{2(dB)} - G_{A(dB)}$$

### 4.3.3 ANTENNA GAIN MEASUREMENT (Method using two identical antennae)

#### ■ Composition :

- 1 GUNN Source : OSG 100
- 1 Ferrite isolator : ISO 100
- 1 PIN Modulator : MOD 100
- 1 Micrometer attenuator : ATM 100
- 1 Calibrated attenuator : ATC 100
- 1 Wavemeter : OND 100
- 2 15 or 20 dB horn antennae : ANC 100/15 or ANC 100/20
- 2 Fast attachment adapters : AFR 100
- 1 Waveguide/coaxial transition : TGN 100
- 1 Measurement sensor : DEN 100 or 1 Thermocouple probe : ST 613
- 4 Waveguide supports : SUP 100
- 1 Measurement indicator : IR 205 or 1 Powermeter : MH 600
- 1 GUNN power supply: CF 204

#### ■ Activation of the bench

The components used for the experiments that follow are essentially BDH R 100 basic bench components.

- Assemble the microwave components and connect the CF 204 GUNN power supply and the IR 205 (or MH 600 and ST 613) measurement indicator as indicated below (fig. 1).
- Power up the GUNN power supply and the measurement indicator.
- Make the measurements.

Identical cases of two horns


fig : 1

The distance R between the two ANC 100 horns must be such that:

$$R \geq 2d^2 / \lambda_0$$

Therefore, it is necessary to measure the frequency on the wavemeter and work out  $\lambda_0$ .

The horns must be accurately aligned on the two planes E and H (fig : 1).

- Make sure there are no obstacles in the horn radiation field and eliminate any spurious reflections using absorbing panels.

■ **Proceed as follows :**

- Set the attenuator ATC 100 to the zero position
- Adjust the ATM 100 attenuator to obtain correct reading on the galvanometer of the IR 205 SWR indicator
- Locate this indication, i.e.  $I_1$
- Eliminate the two horns without modifying the previous setting of ATM 100 and set ATC 100 to 40 dB, the assembly is then as shown in the following figure with connection between the two benches :


fig : 2

- Adjust the ATC 100 calibrated attenuator to the same indication as identified on IR 205 with the previous setup, i.e.  $I_1$ .
- On ATC 100, read the power ratio  $P_r/P_e$  in dB.
- Enter this value into the relation (5) of the theoretical reminders and use the values  $\lambda_0$  and R determined above..

The double of the gain is obtained directly in dB.

Take care to use consistent units: d in m and R in m.

## 4.4. HORN ANTENNA STUDY

### 4.4.1 Introduction

The horn radiation pattern measurement consists in determining :

- the amplitude and shape of the major lobe of the radiation pattern and its width at 3 dB
- the amplitude of the side lobes

For a rectangular horn, therefore using rectilinear polarization, these measurements are made in the two planes E and H, in elevation and azimuth.

### 4.4.2 Theoretical reminders

An isotropic antenna is designed so that a feed point radiates consistently throughout space.

A real antenna preferably transmits in certain directions, making for maximum gain.

A radiation pattern is a representation of the field intensity or of the power transmitted according to the angle of observation at a constant distance.

The radiation pattern of a horn is represented either in plane E (parallel to the smaller side of the waveguide) or in plane H (parallel to the longer side of the waveguide).

The radiation pattern consists of several lobes, and the greater part of the energy is concentrated in the major lobe. Generally, an antenna is designed to minimize the power radiated in the side lobe and the rear lobes to improve directivity

The width of the major lobe at 3 dB is the angle between two points where the power received equals half the maximum value.

### 4.4.3 Radiation pattern measurement

#### ■ Composition :

- 1 Transmission bench :	BEM 100		
- 1 Calibrated antenna :	ATC 100		
- 1 GUNN power supply :	CF 204		
- 1 Rotating platform :	PTM 100	ou	1 rotary joint : JTG 100
- 1 20 dB horn antenna :	ANC 100/20		
- 1 Fastening adapter:	AFR 100		
- 1 Waveguide/coaxial transition :	TGN 100		
- 1 Measurement sensor :	DEN 100	ou	1 Thermocouple probe : ST 613
- 1 Antenna support :	SAN 100		
- 1 Measurement indicator :	IR 205	ou	1 Powermeter : MH 600

#### ■ Activation

This bench consists of an ANC 100/20 dB horn attached to a PTM 100 rotating platform by an SAN 100 antenna support or mounted on a JTG 100 rotary joint.

A TGN 100 waveguide/coaxial transition with a DEN 100 detector, connected to an SWR indicator, is used for making the measurement of the reception level.

Set up the microwave components and connect the CF 204 power supply and the IR 205 SWR indicator as indicated below : (fig. 1)

- Turn the power supply and the SWR indicator
- Make the measurements


fig : 1

■ **Distant field radiation condition:**

Take the following precautions:

- Provide for a minimum distance between the two ANC 100 horns:

$$R \geq 2d^2 / \lambda_0$$

R = distance between transmission antenna and reception antenna

d = largest size of largest antenna

$\lambda_0$  = Wavelength in air

- The horns must radiate in free space, i.e. distant from any metal obstacles, and in particular far from the floor.
- The horn axes must be accurately aligned on the horizontal plane and vertical plane. Check this alignment by moving the reception horn slightly until maximum deflection is obtained on the IR 205 galvanometer.

Proceed as follows :

- Set the ATC 100 attenuator to zero.
- Adjust the IR 205 output level to obtain correct deflection, using the ATC 100 attenuator.
- Locate the attenuation  $A_0$  read and deflection  $D_0$  on 205.
- Turn the receiver horn through an angle  $\alpha_1$ .
- Bring the ATC 100 attenuator to attenuation  $A_1$  until the point  $D_0$  read previously on IR 205 is obtained.
- Turn an angle  $\alpha_2$ , set the ATC 100 attenuator to  $A_2$  to refer back to the  $D_0$  indication on IR 205.

Perform these manipulations for several  $\alpha$ , rotation angle values either side of the horn axis and observe for each angle:

$$\Delta_{1(dB)} = A_0 - A_1$$

$$\Delta_{2(dB)} = A_0 - A_2 \text{ etc...}$$

- Plot the  $\Delta_{dB}$  curve as a function of  $\alpha$  on a diagram using circular coordinates, corresponding to the horn radiation pattern.

#### 4.4.4 Gain Measurement

Refer to paragraph 4.3

## 4.5. DIELECTRIC ANTENNA STUDY

### 4.5.1 Introduction

These longitudinal radiation antennae were imagined and developed by ZINKE and MALLACH. They are used almost only for microwave applications; they have a relatively wide natural bandwidth. Unfortunately, their excitation is relatively unsuitable for broadband setups, limiting their use to special cases.

The purpose of this experiment is to plot the radiation pattern on a dielectric antenna used as a receiving antenna.


### 4.5.2 Theoretical reminders

The concept of a dielectric antenna radiation pattern can be presented by analogy to a radiating wire. Let us consider an energy transmission line composed of elements having a radiation pattern.

$$e(\theta)$$

which we assume, for simplification purposes, to be revolving about the axis of the line.

As an initial approximation, it is accepted that the current along the line is constant and that the line is closed on its characteristic impedance. It is also assumed that the propagation velocity  $v$  over

the line is different from  $c$ , equivalent to considering that the wavelength  $\lambda_g$  is different from  $\lambda_0$ , the wavelength in the air. The line propagation constant is :

$$\alpha = j \frac{2\pi}{\lambda_g} = j \frac{2\pi \lambda_0}{\lambda_0 \lambda_g} = jk \frac{\lambda_0}{\lambda_g} \quad \text{avec} \quad k = Nb \text{ d'onde}$$

and the current at  $M$  is :

$$I = I_0 e^{jk \frac{\lambda_0}{\lambda_g} x}$$

The radiation of the element  $M$  is in advance by  $(kx \cos \theta)$  with respect to the element in  $O$ , and the expression of the field radiated by the line has the form :

$$E = \frac{60 \pi L I_0}{\lambda_0 R} e(\theta) \frac{\sin \pi L \left( \frac{1}{\lambda_g} - \frac{\cos \theta}{\lambda_0} \right)}{\pi L \left( \frac{1}{\lambda_g} - \frac{\cos \theta}{\lambda_0} \right)}$$

where  $L$  is the antenna length

In this expression the fraction is in the form  $\frac{\sin u}{u}$  and this function represents a maximum for  $u = 0$ .

The maximum radiation will therefore occur for

$$\frac{1}{\lambda_g} = \frac{\cos \theta}{\lambda_o}$$

Thus:

$$\cos \theta = \frac{\lambda_o}{\lambda_g}$$

Therefore, if  $\lambda_g$  is the order of  $\lambda_o$ , there will be maximum radiation in the direction of the axis. It will be seen that if the line is a simple metal wire, then :

$$e(\theta) = \sin \theta$$

giving zero in the direction of the wire.

### 4.5.3 Radiation pattern measurement

#### ■ Composition

- 1 Transmission bench :	BEM 100	
- 1 Turntable :	PTM 100	
- 1 Antenna support :	SAN 100	
- 1 Dielectric antenna :	AND 100	
- 1 Measurement sensor :	DEN 100	or 1 Thermocouple probe : ST 613
- 1 Waveguide/Coaxial transition :	TGN 100	
- 1 Micrometer short-circuit :	CCM 100	
- 1 Measurement indicator :	IR 205	or 1 Powermeter : MH 600
- 1 GUNN power supply :	CF 204	

#### ■ Activation of the reception dielectric antenna bench

This bench consists of one or several AND 100 dielectric antennae terminating at one end by a CCM100 micrometer variable short-circuit and at the other by a TGN100 waveguide/coaxial transition equipped with a DEN100 detector which, connected to an IR205 SWR indicator, allows reception level measurement to be made.(fig 2)

These AND 100 antennae can be mounted in arrays of two, three or four, as indicated in the drawing below (Fig. 1), to obtain particularly pronounced directivity; they are then powered in parallel and in phase agreement.


fig : 1

- Assemble the microwave components forming the reception bench and connect up the IR 205 SWR indicator as indicated below (fig.2).
- Mount this reception system on the PTM 100 rotating platform.
- Power up the CF 204 GUNN power supply and the IR 205 SWR indicator.
- Adjust the CCM 100 short-circuit micrometer to obtain maximum deflection on the IR 205 indicator.
- Proceed with the measurements.


fig : 2

■ **Radiation pattern measurement**

Proceed as described in paragraph 4.7.3 of the parabolic antenna study.

**4.5.4 Gain measurement**

Proceed as described in paragraph 4.3.

## 4.6. STUDY OF SLOT ANTENNA

### 4.6.1 Introduction

A slot antenna consists of a straight waveguide section whose wall corresponding to the large side has several slots arranged so that they radiate.

Slots blanking pieces are used for varying the width of the radiation lobe according to the number of active slots.

The purpose of this experiment is to plot the radiation patterns corresponding to various configurations: active slots numbers, two, five, then ten, included in the antenna.

### 4.6.2 Theoretical reminders

Slot-shaped openings in the wall of a guide radiate if they intersect the current lines.

The field radiated by this type of slot is similar to that of a thin dipole having the same dimensions (Babinet Principle) - (paragraph 4.6.3). Therefore, the electrical field is expressed similarly.

#### ■ Simple slot radiation pattern

With a slot placed at the origin of a system of axes (o, x, y, z), aligned along o x, on the o z axis and in the y o z plane, the field will be perpendicular to o z (fig : 1).


fig : 1

If we move away from this angle by an angle  $\theta$ , the electrical field amplitude will be:

$$E = EM \frac{\cos \left[ \frac{\pi}{2} \sin \theta \right]}{\cos \theta}$$

In the x o z plane, the pattern is identical to that of the dipole (fig. 2).


fig : 2

It is revolving with respect to the o x axis. We will only consider the half-space  $z > 0$  in the following, and the radiation will be observed in the plane perpendicular to the slot.

■ **Equidistant slot radiation pattern**

Let us consider  $n$  identical radiating slots, very thin and regularly spaced at an interval  $d$  (fig. 3).


fig : 3

These slots are illuminated with an angle of incidence  $\theta$ , and together form  $n$  sources,  $A_0, A_1, A_{n-1}$ , regularly spaced at a distance  $d$  (fig. 4).


fig : 4

The step difference is  $(d \sin \theta)$ , and the phase-shift between the waves  $B_0 A_0$  and  $B_1 A_1$ , captured respectively by  $A_0$  and  $A_1$  is:

$$\phi = 2 \pi \sin \theta$$

Thus  $E_0 = E e^{j\omega t}$ , field in  $A_0$

The field in  $A_1$  is written:  $E e [j (\omega t - \phi)]$

The field in  $A_{n-1}$  is written:  $E e \{ j [(\omega t - (n-1) \phi)] \}$

If there is no phase-shift between feeds  $A_0, A_1 \dots, A_{n-1}$ , replacing  $\phi$  by its value, and considering the field modulus and a small  $\theta$ , the resultant field is given by :

$$|E_r| = |E| \frac{\sin \left[ \frac{n \pi d}{\lambda} \theta \right]}{\sin \left[ \frac{\pi d}{\lambda} \theta \right]}$$

If  $(2\pi d \theta / \lambda)$  is small, the sine can be identified by its argument and :

$$|E_r| = |E| \frac{\frac{n \pi d \theta}{\lambda}}{\frac{\pi d \theta}{\lambda}} \approx n |E|$$

For  $\theta \neq 0$ , the total field is canceled for values of  $\theta_n$ , such that:

$$\frac{n \pi d}{\lambda} \theta_n = k \pi \quad \theta \neq 0$$

In particular, for  $k = 1$ , we have:

$$\frac{n \pi d}{\lambda} \theta_1 = \pi$$

$\theta_1$  is the first value for which  $E_r$  is zero, determining the width of the major lobe.

Assuming that :

$$\theta_1 = \lambda / nd$$


It can be seen that  $\theta_1$  is increasingly small in proportion to the high number  $n$  of feeds used. If  $n$  is high, and in particular for  $\lambda \ll n\theta_1$

$$\sin \left[ \frac{\pi d}{\lambda} \theta \right] \approx \frac{\pi d}{\lambda} \theta$$

$$E(u) = n E \frac{\sin u}{u}$$

with:  $u = n (\pi d / \lambda) \theta$


We obtain the conventional curve (fig. 5).

Note that the slots used are of the "shunt" type, i.e. parallel to the longitudinal axis of the waveguide. These slots are less sensitive to the frequency than the transversal slots.

### 4.6.3 Babinet principle

Let us consider two planes, one of which is plane A, in which screens are arranged, and the other plane B, back from A with respect to any light source S, and in which we observe the field.

In this case (1), we interpose an absorbing screen in A, so that the field on B depends on :

$$F_e = f_1(x, y, z)$$

In the case (2), we interpose a complementary screen in A, so that the field on B is now :

$$F'_e = f_2(x, y, z)$$

The **BABINET** Principle specifies that if the field without a screen (case 3) is :

$$F_o = f_0(x, y, z)$$


we still have :

$$f_1(x, y, z) + f_2(x, y, z) = f_0(x, y, z)$$

This is obvious for geometrical optics, as shown by simple examination of the figure, but the principle also applies in the case of diffraction.

However, here, no allowance is made for the vector property of electromagnetic radiation, so that the BABINET Principle had to be generalized by BOOKER, who demonstrated that if the screen considered is a perfect electrical conductor, its complementary screen must be a perfect magnetic conductor.

A conductor like this does not exist, but the same effect can be obtained using two perfect electrical conductors (e.g. copper screens) while switching the electrical and magnetic fields.


For instance, let there be a feed S radiating a horizontal electric field and a metal screen M producing in the Q plane a field  $E_1$  (case 1). The complementary screen will be an aperture having the same size as the first screen, pierced in a metallic plane (case 2), but the polarization will have to be turned through  $90^\circ$  to replace the electric fields by magnetic fields.

We can also use an aperture having the same size again, but arranged normally with respect to the original screen (case 3) so that polarization is preserved.

On the basis of these considerations, BOOKER also demonstrated that the field radiated by an aperture 0 (case 2 above) was the same as radiated by a thin dipole having the same size as the aperture (screen M, case 1 above), excited at its center and operating with polarization at right angles to the aperture.

Accordingly, the fields radiated by the two devices in the figure below are equal in amplitude, but the electric field polarizations are perpendicular :

#### 4.6.4 Radiation pattern measurement

##### ■ Composition :

- |   |         |    |                              |
|---|---------|----|------------------------------|
| - 1 Transmission bench :                | BEM 100 |    |                              |
| - 1 Rotating platform :                 | PTM 100 |    |                              |
| - 1 Slot antenna :                      | ANF 100 |    |                              |
| - 1 Micrometer variable short-circuit : | CCM 100 |    |                              |
| - 1 Waveguide/coaxial transition :      | TGN 100 |    |                              |
| - 1 Measurement sensor :                | DEN 100 | or | 1 Thermocouple probe: ST 613 |
| - 1 Antenna support :                   | SAN 100 |    |                              |
| - 1 Measurement indicator :             | IR 205  | or | 1 Powermeter: MH 600         |
| - 1 GUNN power supply :                 | CF 204  |    |                              |

##### ■ Activation of the "slot" antenna reception bench

This bench consists of an ANF100 slot antenna, closed at one end by a CCM 100 variable short-circuit, and at the other by a TGN 100 waveguide/coaxial transition equipped with a DEN 100 detector, which, connected to an IR 205 SWR indicator, allows the reception level to be measured.

The distance to separate transmission from reception must be greater than 9 meters (see paragraph 3.1).

- Assemble the microwave components forming the reception bench with connection of the SWR indicator as indicated below. (fig 1)
- Attach this reception assembly to the PTM 100 rotating platform
- Power up the CF 204 GUNN power supply and the IR 205 SWR indicator.
- Adjust the CCM 100 short-circuit micrometer to obtain maximum deflection on the IR 205 indicator.
- Proceed with the measurements.


fig.1

### ■ Radiation pattern

This manipulation consists in measuring the general form of the major lobe radiation pattern between the first two side lobes "LS<sub>d2</sub>" et LS<sub>g2</sub>".

The measurement of the radiation levels received for sufficiently small angular deviations is made to obtain good radiation pattern trace definition.

**Note:** To make this measurement, it is advisable to use the "extended" scale of the IR 205 level indicator.

This provides greater accuracy, more particularly necessary near the axis of symmetry, because variations are subsequently far greater.

- Place the zero graduation, the rotating platform origin, on the transmission feed axis of symmetry and the slot antenna axis.
- Check that the maximum deflection of the identifier corresponds to the rotating platform origin "Zero" by moving the antenna either side of the axis of symmetry.
- Adjust the rotating platform "Zero" index for this maximum, or adjust the alignment of the reception bench to obtain this condition.
- Set the level indicator gain.

### ■ Pattern with two slots

Blank off all the slots except for the two at the center of the antenna.

Make measurements on the radiation pattern:

- side lobe width at 3 dB,
- side lobe level,
- radiation pattern.
- major lobe width at 3 dB
- turn the rotating platform to the right to obtain 3 dB attenuation with respect to the reference 0 dB on l'IR 205.
- make a note of angle  $\theta_{g1}$  on the rotating platform facing the mark,
- return to zero.
- turn the rotating platform to the left to obtain 3 dB attenuation and make a note of angle  $\theta_{d1}$  of the rotating platform facing the mark;
- return to zero.
- the angle  $\theta_{d1}$  must be identical to the angle  $\theta_{g1}$
- the major lobe width at -3 dB equals:  $\theta_{d1} + \theta_{g1}$

### ■ Side lobe level

Proceed in the same way as for measuring the major lobe at - 3 dB and continue to turn the antenna support rotating platform.

The level of the received signal passes through a minimum, then through a maximum; the first maximum corresponds to the level of the first side lobe; by continuing the rotation, pass successively through a minimum, then a maximum to obtain the following side lobe levels.

Depending on the direction of rotation of the rotating platform (to the right or left), we locate the first right side lobe "LS<sub>d2</sub>" or left side lobe "LS<sub>g2</sub>", then the following "LS<sub>d3</sub>" or "LS<sub>g3</sub>", etc..

Make a note of the angular deviations of each of the side lobes :

- " $\theta_{d2}$ " corresponding to "LS<sub>d2</sub>" and " $\theta_{g2}$ " corresponding to "LS<sub>g2</sub>"
- " $\theta_{di}$ " corresponding to "LS<sub>di</sub>" and " $\theta_{gi}$ " corresponding to "LS<sub>gi</sub>"

**Note** The received signal level for each of these side lobes is obtained by varying the sum of the attenuation value obtained on the IR 205 SWR indicator (keypad + deviation).

The results :  $X_{(dB)} = f(\theta)$

can be represented in two different forms :

■ In Cartesian coordinates


fig : 2

■ In polar coordinates



fig : 3

■ Pattern with five slots

Blank off the slots except for 5 of them either side of the antenna center. Make the measurements in the same way as before.


Slot antenna :  
5 slots uncovered  
toward short-circuit

fig : 4

■ **Pattern with all slots**

Make the measurements in the same way as for the pattern with 5 slots (previous paragraph)

Slot antenna:  
all slots

fig : 5

**4.6.5 Gain measurement**

Proceed as described in paragraph 4.3.

## 4.7. PARABOLIC ANTENNA STUDY

### 4.7.1 Introduction

The parabolic antenna, using a parabolic reflector, is widely used for microwave applications. Like a converging lens for optical purposes, it transforms the spherical waves from a source (known as primary) at its focal point into a beam of parallel plane waves, to obtain very high directivity. It is even stronger if the primary feed is the horn because of the relative decrease of the side lobe intensity.

The purpose of this experiment is to plot the radiation pattern of a parabolic antenna used for reception and to locate the higher order lobes.

### 4.7.2 Theoretical reminders

#### ■ Paraboloid of revolution properties

The paraboloid is a reflecting element whose characteristics are widely used for microwave energy radiation.

The purpose of this paragraph is to reveal the influence of its characteristics and the positioning of the feed on antenna radiation.

An individual antenna placed at the focal point of a paraboloid produces a beam of parallel rays by reflection.


Although the diffraction phenomenon predominates in the microwave field, this main property of the paraboloid demonstrates that in the "P" plane, at right angles to its axis, the wave has an equiphase surface area according to the figure below and the accompanying relation.

The equation of the paraboloid is:  $y^2 = 4 f x$ ,  
f being the focal length of the parabola

$$FM + MM' = FN + NN' = 2f$$

$$f = \frac{D}{4} \left[ \frac{1 + \cos \phi_0}{\sin \phi_0} \right]$$

$$D = 4f \frac{\sin \phi_0}{1 + \cos \phi_0} = 4f \operatorname{tg} \frac{\phi_0}{2}$$

For microwave applications, the energy is shared by diffraction in the form of a major lobe. This is accompanied by side lobes with levels varying with the distribution on the reflector and the attenuation on the edge of the latter.

#### ■ Expression of parabolic antenna gain

At a certain distance, an antenna can be considered as an individual source whose power radiated by a solid unit of angle depends on the direction. If the latter is defined by the angles  $\theta$  and  $\phi$ , the antenna gain will be given in the following relation :

$$G_{(\theta, \phi)} = \frac{P_{(\theta, \phi)}}{P_0 / 4\pi}$$

where :  $\theta$  elevation angle (horizontal)  
 $\phi$  azimuth angle (vertical)  
 $P_{(\theta, \phi)}$  is the power radiated in this direction  
 $P_0$  is the total radiated power  
 $P_0 / 4\pi$  is the power per solid unit of angle

In general, the gain of an antenna represents the concentration of the power radiated in a privileged direction with respect to the power that a radiating isotropic feed  $P_0$  would radiate in this direction.

In the particular case of a paraboloid of revolution, the maximum gain direction coincides with the direction of its major axis :

$$G_{(\theta,\phi)} = \frac{P_{(\theta,\phi)}}{P_o / 4\pi}$$

However, for a parabolic reflector,  $G_{(0,0)}$  is given by:

$$G_{(0,0)} = \left[ \frac{\pi D}{\lambda} \right]^2 \left[ \int_0^{\phi_0} \sqrt{\gamma_{(0)}} \operatorname{tg} \frac{\phi}{2} d\phi \right]^2 \operatorname{cotg}^2 \left( \frac{\phi_0}{2} \right)$$

where  $\gamma_{(0)}$  is the gain of the primary source in the direction  $\phi$ .

The maximum theoretical gain of a uniformly illuminated circular aperture is is:

$$G = \frac{4\pi \Sigma}{\lambda} = \frac{4\pi}{\lambda} \left[ \pi \frac{D^2}{4} \right] = \frac{\pi^2 D^2}{\lambda^2}$$

$\Sigma$  is the equivalent surface area of the radiating aperture.

We can therefore deduce the theoretical gain factor F

$$F = \frac{G_{(\theta,\phi)}}{(\pi D / \lambda)^2} = \left[ \int_0^{\phi_0} \sqrt{\gamma(\theta)} \operatorname{tg} \frac{\phi}{2} d\phi \right] \operatorname{cotg}^2 \left( \frac{\phi_0}{2} \right)$$

This relation does not take into consideration losses due to the energy that is not intercepted by the reflector (Spill-over).

Therefore, to obtain the real gain factor, we must multiply  $F_0$  by the illumination efficiency « D » :

$$\eta = P_i / P_p$$

$P_i$  the power intercepted by the paraboloid


$P_p$  the primary emitted power

### ■ Radiation pattern

The radiation pattern is a curve representative of the relative gain as a function of the angle  $\theta$ . This measurement is made in a horizontal plane containing the axis of symmetry of the paraboloid.

The relative gain is 1 on the axis of symmetry and less than 1 everywhere else.

The relative gain is generally evaluated in dB; the value of 0 dB, the reference value, is taken on the axis of symmetry.

■ **Distant field radiation condition:**

Comply with the minimum distance between the two ANC 100 horns:

$$R \geq 2d^2 / \lambda_0$$

R = distance between transmission antenna and reception antenna

d = largest size of the largest antenna

$\lambda_0$  = Wavelength in air

Any assembly to be installed for testing comprises a transmission bench and a reception bench. The distance L or R separating them is linked with the size of the antennas and the  $\lambda$  wavelength transmitted by the transmission source.

The phase difference between the wave affecting the edge of the receiving antenna and the wave in the axis must be at most equal to  $\lambda/16$ . This remains within the distant field radiation conditions.

$$\Delta R = R' - R$$

$$d^2 / 4 + R = (R + \Delta R)^2$$

Hence:

$$R^2 + 2R\Delta R + (\Delta R)^2 = R^2 + d^2 / 4$$

Leaving out the second order term:


$$R = d^2 / 8\Delta R$$

Well:

$$\Delta R \leq \lambda / 16$$

Hence

$$R \geq 2d^2 / \lambda$$

### 4.7.3 Radiation pattern measurement

■ **Composition :**

- 1 Transmission bench: BEM 100
- 1 Rotating platform: PTM 100
- 1 Parabolic reflector: ANP 100
- 1 10 dB horn antenna: ANC 100/10
- 1 Fastening adapter: AFR 100
- 1 Waveguide/Coaxial transition: TGN 100
- 1 Measurement sensor: DEN 100 or 1 Thermocouple probe : ST 613
- 1 GUNN power supply: CF 204
- 1 Measurement indicator: IR 205 or 1 Powermeter : MH 600

■ **Activating the parabolic reception bench**

This bench consists of an ANP100 parabolic reflector having diameter 350 mm. The reception circuit consists of a horn (10 dB Gain), ANC 100/10, placed at its focal point on its main axis. A TGN 100 waveguide/coaxial transition provided with a DEN 100 detector, connected to an IR205 SWR indicator, allows the reception level measurement to be made.

The supplied ANP100 paraboloid, having diameter d = 350 mm, results in wavelength  $\lambda_0$  being 30 mm approximately, with a minimum distance R between the aperture of the transmission horn and the reception antenna of :

$$R \geq 2d^2 / \lambda_0 \approx 9m$$

The distance that has to separate transmission from reception must be greater than 9 meters,

Assemble the microwave components forming the reception bench and connect up the IR 205 SWR indicator as shown below (fig. 2).

- Attach this reception assembly to the PTM 100 rotating platform
- Power up the CF 204 power supply and IR 205 SWR indicator.
- Make the necessary measurements.


Fig. 2.

■ **Measurement of radiation pattern**

Set the rotating platform origin graduation on the axis of symmetry of the transmission source and the axis of the paraboloid (fig. 3).


fig. 3

**Caution:**

Check that the maximum deviation of the IR 205 indicator corresponds to the rotating platform origin "zero" by moving the paraboloid on either side of the axis of symmetry. Adjust the index of the rotating platform to "zero" for this maximum, or adjust the alignment of the reception bench to obtain this condition. Adjust the level indicator gain.

■ **Width of the major lobe at 3 dB**

- Turn the rotating platform to the right to obtain 3 dB attenuation.
- Make a note of the angle  $\theta_{d1}$  of the rotating platform facing the mark.
- Return to zero.
- Turn the rotating platform to the left until 3 dB attenuation is obtained.
- Make a note of the angle  $\theta_{g1}$  of the rotating platform facing the mark.
- Return to zero.

The angle  $\theta_{d1}$  must be identical to the angle  $\theta_{g1}$ .

The width of the major lobe at -3 dB equals  $\theta_{d1} + \theta_{g1}$

■ **Side lobe level**

Proceed in the same way as for major lobe width measurement at – 3 dB and continue to turn the antenna support rotating platform.

The level of the received signal passes through a minimum, then increases again to pass through a maximum. This maximum corresponds to the side lobe level.

Depending on whether the rotating platform turns to the right or to the left, locate the first right side lobe  $LS_{d2}$  or the first left side lobe  $LS_{g2}$ .

- Make a note of the angular deviation between each of the side lobes  $\theta_{d2}$  and  $\theta_{g2}$ .
- Make a note of the signal level received for each of these side lobes by working out the sum of the value of attenuations obtained on IR 205 (those of the attenuation switch and the value read on the level indicator galvanometer).

**Note:** In this way, it is possible to locate the following side lobe levels:  $LS_{d2}$ ,  $LS_{g2}$ ,  $LS_{d3}$ ,  $LS_{g3}$ , etc... while continuing to turn the rotating platform and making measurements in the same way.

■ **Radiation pattern**

The process involves measuring the general form of the major lobe radiation pattern between the two first side lobes  $LS_{d2}$  and  $LS_{g2}$ .

In this case, it is necessary to measure the radiation levels received for sufficiently small angular deviations to permit high definition of the radiation pattern trace.

To make this measurement, it is advisable to use the "extended" scale of the level indicator. It provides greater precision, which is more particularly necessary near the axis of symmetry because the variations are subsequently far greater.

A pattern as indicated in the figure (Fig. 3) is obtained by entering the measured points onto a circular pattern.


fig : 4

**4.7.4 Gain measurement**

Proceed as described in paragraph 4.3.

## 4.8. PASSIVE RADAR RESPONDER STUDY

### 4.8.1 Introduction

In 1944, during his course on optics at Brown University, Professor R.K. Lüneburg described the general law enabling a spherical lens with continuous index variation to have a given image correspond to a given source. The formula is:

$$\varepsilon = 2 - (r / r_0)^2$$


where  $\varepsilon$  = permittivity of the environment

$r_0$  = radius of the sphere

$r$  = current radius of the sphere (variable value included between 0 and  $r_0$ )

The Lüneburg lens, a direct application of this law, led to the Lüneburg reflector used as passive radar responder.

Today, **LUN'TECH** France manufactures and markets these dielectric lenses whose electromagnetic and mechanical characteristics correspond to the needs of military applications, such as the radar identification of missiles and targets, and develops civil applications like the beaconing of ports and airports, vehicle, ship signaling, buoys, obstacles, etc..

### 4.8.2 Theoretical reminders

Lüneburg lenses are currently used as reflectors in passive radar responders. They increase substantially the Radar Equivalent Surface, "R.E.S.", (several times greater than that of a metal sphere having the same diameter) for an object containing this type of reflector. There are monostatic reflectors, reflecting all the energy received into a cone having a small angle of aperture, and bistatic reflectors that reflect energy into a wider cone (on the order of  $\pm 10^\circ$ ).

The Lüneburg reflector consists of a Lüneburg lens at the focal point of which there is a spherical metal surface (fig 1).


fig : 1

#### ■ Radar Equivalent Surface (R.E.S.)

##### Definition :

The Radar Equivalent Surface translates the aptitude of the lens to return radar waves in a given direction..

A reflector is considered as an antenna with a surface area  $S$  and a gain  $G$ . Its RES is then defined by the product of these two magnitudes in a given direction [9] :

$$\sigma = S.G \quad (\text{in m}^2).$$

When antennae are large compared to the wavelength, their surface tends to bend into their main structure.

The relation between the gain and surface of an antenna is as follows:

$$G = \eta \cdot \left[ \frac{4 \cdot \pi S}{\lambda^2} \right] \text{ with } 0 < \eta < 1$$

For a reflector, the RSE equals the product of the equivalent gain G by the effective surface of the isotropic antenna.

$$\sigma = G^2 \cdot \frac{\lambda^2}{4 \cdot \pi}$$

i.e :

$$\sigma = \eta^2 \cdot \frac{4 \cdot \pi \cdot S^2}{\lambda^2}$$

thus, for a revolving reflector having diameter D :

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



As an example, here is the RSE of several metallic reflectors in the X band (at a frequency of 9,4 GHz).

Shape	Dimensions (b=R=20 cm)	Theoretical RSE <sup>1</sup> $4 \pi S^2 / \lambda^2$	Notes
Metallic trihedron		$4 \pi ((3 \sqrt{3}b^2) / 8)^2 / \lambda^2 = 8.28 \text{ m}^2$	Average radar surface Response angle greater than for disk
Metallic disk		$4 \pi (\pi R^2)^2 / \lambda^2 = 195.0 \text{ m}^2$	Large radar surface identical to that of Lüneberg reflector Very fine response angle
Metal sphere		$\pi R^2 = 0.12 \text{ m}^2$	Small radar surface Isotropic
Lüneburg spherical dielectric reflector		$4 \pi (\pi R^2)^2 / \lambda^2 = 7.75 \cdot D^4 / \lambda^2 = 195.0 \text{ m}^2$	Large radar surface Wide response angle

<sup>1</sup> The RSE of a sphere is independent of a frequency if the wavelength is small compared to its size, i.e.  $\lambda \ll R$ .

■ Link check-up


fig : 2

- A : Amplifier.
- Pr : received power.
- Pe : transmitted power.
- Ge : transmission horn gain = 15 dB.
- Gr : reception horn gain = 15 dB.
- Pespace : losses in space.
- Pcables : losses in cables.
- D : Distance between antennae and dish.

The balance of the link in decibels works out as follows [8] :

$$P_r = P_e + G_e + G_r + P_{espace} + \text{contribution du disque} - P_{cables}$$

$$P_r = P_e + G_e + G_r + 10 \cdot \log \left\{ \left( \frac{\lambda}{4 \cdot \pi \cdot D} \right)^2 \right\} + 10 \cdot \log \left( \frac{\sigma}{4 \cdot \pi \cdot D^2} \right) - P_{cables}$$

However

$$\sigma = SG = \frac{4 \pi S^2}{\lambda^2}$$

Therefore

$$\frac{G}{4 \pi D^2} = \frac{S^2}{\lambda^2 D^2}$$

Hence,

$$P_r = P_e + G_e + G_r - P_{cables} + 10 \cdot \log \left\{ \left( \frac{\lambda}{4 \cdot \pi \cdot D} \right)^2 \right\} + 10 \cdot \log \left( \frac{S^2}{\lambda^2 D^2} \right)$$

■ Numerical application

For  $f = 9.375$  GHz and a disk having a diameter of 127 mm (radius = 63.5 mm) transmission horn gain ( $G_e$ ) : 15 dB. reception horn gain ( $G_r$ ) : 15 dB.

transmitted power ( $P_e$ ) : 13 dB. losses in cables ( $P_{cables}$ ) : 5 dB.

Distance between antennae and dish (D) : 13 m.

$$P_r = 38 + 10 \cdot \log \left\{ \left( \frac{3,2 \cdot 10^{-2}}{4 \cdot \pi \cdot 13} \right)^2 \right\} + 10 \cdot \log \left( \frac{(\pi \cdot (6,35 \cdot 10^{-2}))^2}{(3,2 \cdot 10^{-2})^2 \cdot 13^2} \right)$$

$$P_r = 38 - 74,16 - 20,38$$

$$P_r = -56,54 \text{ dBm}$$

■ **Relations used during measurements**

The measurements are made comparatively with a metallic disk.  
The linear link breakdown for a disk is following, without allowing for losses in the cables:

$$P_{Rd} = P_e G_e G_r \left( \frac{\lambda}{4\pi D} \right)^2 \frac{\sigma_d}{4\pi D^2} = \frac{P_e G_e G_r \lambda^2}{(4\pi)^3 D^4} \sigma_d$$

For the reflector, the relation is identical, but  $P_{Rd}$  is replaced by  $P_{Rr}$  and  $\sigma_d$  by  $\sigma_r$ .

These relations show that the power ratio equals the radar surface ratio:

$$\frac{P_{Rr}}{P_{Rd}} = \frac{\sigma_r}{\sigma_d}$$

During change over to decibels, this relation gives:

$$10 \log(P_{Rr}) - 10 \log(P_{Rd}) = 10 \log\left(\frac{\sigma_r}{\sigma_d}\right)$$

Therefore

$$\frac{P_{Rr}|_{dB} - P_{Rd}|_{dB}}{10} = \log\left(\frac{\sigma_r}{\sigma_d}\right)$$

Giving the following relation used during measurements:

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- with :  $\lambda$  the wavelength (in m).
- $SER_{référence}$  : the theoretical RSE of the reference disk :
- $G_{référence}$  : the level with the measured reference disk
- $G$  : the level measured with the investigated reflector.

### 4.8.3 Measurement of RSE (Equivalent Radar Surface)

■ **Composition :**

- 1 Transmission bench : BEM 100
- 1 15 dB horn antenna : ANC 100/15
- 1 Fastening adapter : AFR 100
- 1 Waveguide/coaxial transition : TGN 100
- 1 Measurement sensor : DEN 100 or 1 Thermocouple probe: ST 613
- 1 Rotating Platform : PTM 100
- 1 Waveguide support : SUP 100
- 1 Antenna support : SAN 100
- 1 Passive radar responder : RRL 100
- 1 Reflector disk : DR 100
- 1 Measurement indicator : IR 205 or 1 Powermeter: MH 600
- 1 GUNN power supply : CF 204

■ **Activation of the bench**

- This bench consists of a transmission assembly, an RRL 100 Lüneburg reflector, a DR 100 reflector disk and a reception system including: an ANC 100/15 horn antenna with a TGN 100 transition and a DEN 100 detector for detecting the reflected wave.
- Set up the microwave components while connecting the CF 204 GUNN power supply and the IR 205 SWR indicator as shown below (fig.3).
  - Power up the CF 204 GUNN power supply and the IR 205 SWR indicator.
  - Proceed with the manipulations.


fig : 3

### ■ Radar equivalent surface (RSE) measurement

- Set the GUNN oscillator to approximately 9 GHz.
- Place the disk on the rotating platform, and on the SWR indicator, locate the maximum deviation in elevation and azimuth to determine the optimal position. The disk surface must be absolutely at right angles to the direction of the reception horn so that all the energy is captured.
- The proceed in the same way as for the measurement of an antenna radiation pattern while varying the elevation either side of the maximum in  $1^\circ$  steps, then plot the curve (1) (fig 4).
- Replace the disk by the Lüneburg lens (the sphere), positioning the "metallic" side, (identified by a label), facing the transmission horn. In the same way as previously, plot the curve
- Turn the Lüneburg lens through  $180^\circ$  (normal use of Lüneburg lens) and plot the curve (2) (Fig. 4) as indicated previously.
- Compare the curves obtain and calculate the RSEv.

**NOTE:** It is particularly important to preserve the "transmitter, receiver" distance constant from one manipulation to the other in order to remain in phase.

Reflector and disk radar surface :

1 = Disk

2 = Lüneburg lens


fig 4

## 4.9. STUDY OF DOPPLER EFFECT

### 4.9.1 Introduction

The purpose of this experiment is to qualitatively, then quantitatively study the DOPPLER effect for microwave applications and the relation between the DOPPLER frequency and the linear speed of displacement of the target, by various means of measurement.

The experiment can be run using a real mobile unit or a SDD 100 displacement simulator, which has the advantage of supplying uninterrupted simulation of the speed of displacement while remaining static.

#### NOTE

When the DOPPLER simulator bench is operating, spurious reflection may interfere with the experiment, and an oscilloscope will detect the unwanted DOPPLER frequencies, recognized by their lack of stability.

Therefore, it is advisable to put the DOPPLER experiment bench in an area suffering from little interference from moving objects and, if possible, having few reflecting surface areas.

Use ABS 100 absorbing panels to decrease the effect of the spurious reflections.

Avoid, for instance, people moving around the laboratory during the experiment and remember the effect of spurious echoes that may be used by lighted fluorescent tubes.

### 4.9.2 Theoretical reminders

#### ■ Definition

The DOPPLER effect is a frequency drift induced by the displacement of a mobile body on the echo of a signal transmitted by a fixed source.

#### ■ Applications

The best known are:

- in the acoustic field:
  - \* The change of the whistling pitch of a train or the siren of a police vehicle moving at high speed.
- In the microwave field:
  - \* Radar, used by the police for speed checks.
  - \* Radar used by alarm systems to ensure volumetric protection.

#### ■ Theoretical reminders

Let there be a source transmitting a sinusoidal wave having frequency "f" and pulsation " $\omega$ ", toward a moving target "M"..

The source and the movable target are far between a distance D.


fig : 1

The time interval  $\delta t$ , needed for the transmitted wave to cover the source-target distance at the speed of light  $c$  is :

$$\delta t = D / c$$

With the target in displacement, distance  $D$  varies according to time.

If the wave transmitted by the source is in the form :

$$f(t) = \cos \omega t$$

The signal received by the target,  $\delta t$  later with

$$\delta t = D / c$$

will correspond to signal  $f'(t)$  transmitted by the source :

$$f'(t) = \cos \omega \left( t - \frac{D}{c} \right)$$

i.e.

$$f'(t) = \cos \left( \omega t - \omega \frac{D}{c} \right)$$

Part of this signal, the echo, will be reflected toward the receiver;

During the return journey, the wave will go through a second phase-shift identical to the forward path phase-shift because the journey times will be identical.

The return signal, available at the receiver, is therefore in the form :

$$f''(t) = \cos \left( \omega t - 2\omega \frac{D}{c} \right)$$

The target speed of displacement " $v$ " is constant, and very small compared to  $c$  ; so that if  $D_0$  represents the source-target distance at moment " $t_0$ ", the distance covered between  $t$  and  $t_0$  , can be expressed as follows :

$$D = D_0 \pm v(t - t_0)$$

The echo time function can be written at moment  $t$ , by bringing in the speed and direction of displacement of the target with respect to the source, in the following form :

$$f''(t) = \cos \left( \omega t - 2\omega \frac{D_0 \pm v(t - t_0)}{c} \right)$$

thus

$$f''(t) = \cos \left( \omega t - 2\omega \frac{D_0}{c} \pm 2\omega v \frac{t}{c} \pm 2\omega v \frac{t_0}{c} \right)$$

Note that in the previous relation, the expression  $2\omega v / c$  represents the phase rotation " $\omega_d$ " brought in by the displacement of the target :

so we are able to write

$$f''(t) = \cos \left\{ (\omega \pm \omega_d)t - \left( 2\omega \frac{D_0}{c} \right) \pm \omega_d t_0 \right\}$$

with

$$\omega_d = 2\omega v / c \quad (1)$$

and

$$f_d = \omega_d / 2\pi \quad . (2)$$

The frequency " $f_d$ " is the DOPPLER FREQUENCY induced by the displacement of the target;

(1) and (2) allows us to write the relation

$$f_d = 2fv / c$$

which links:

The DOPPLER frequency	$f_d$
with the transmission frequency	$f$
with the speed of displacement	$v$
with the speed of light	$c$

### 4.9.3 Speed of displacement measurement

#### ■ Composition :

- |                                   |                  |
|-----------------------------------|------------------|
| - 1 GUNN diode oscillator :       | OSG 100          |
| - 1 Ferrite isolator :            | ISO 100          |
| - 1 Wavemeter :                   | OND 100          |
| - 1 Waveguide parallel detector : | DEG 100          |
| - 1 20 dB horn antenna :          | ANC 100/20       |
| - 2 Waveguide support :           | SUP 100          |
| - 1 Oscilloscope :                | MTX 3252         |
| - 1 GUNN power supply :           | CF 204           |
| - 1 Displacement simulator :      | SDD 100 (Option) |
| - 1 Fastening adapter :           | AFR 100          |

#### ■ Activating the bench

The bench consists of a GUNN source having a sinusoidal frequency OSG 100 for transmission, combined with a DEG 100 waveguide parallel detector for detecting the DOPPLER frequency which will be displayed on an oscilloscope. An ANC 100/20 20 dB horn is used as transmitting and receiving antenna. An ISO 100 ferrite isolator between the oscillator and the detector prevents load variations caused by the latter from affecting the oscillator.

An OND 100 wavemeter placed after the source will be used for measuring the frequency (fig 2).

An SDD 100 displacement simulator (optional) will be placed on the transmitted wave path.

- Assemble the microwave components and connect up the CF 204 GUNN power supply and the oscilloscope as indicated below.
- Turn on the power supply and the oscilloscope.
- Make the measurements.


fig : 2

#### ■ Manipulations

- Adjust the GUNN oscillator to 9 GHz.
- Adjust the vertical sensitive of the oscilloscope to 50mV/cm with continuous coupling to measure the voltage supplied to the detector.
- Adjust the frequency with the OND 100 wavemeter at 9 GHz and adjust the GUNN oscillator to obtain absorption on the oscilloscope screen when, at this exact moment, the frequency is 9 GHz.
- Move the wavemeter by at least 50 MHz.
- Start up the SDD 100 displacement simulator.

- Increase the sensitivity of the oscilloscope (AC coupling) until it is possible to study the sinusoidal trace of the DOPPLER frequency. By choosing the right scan speed, measure the period of the signal and work out the corresponding DOPPLER frequency. On the basis of this information, calculate the speed of displacement of the mobile unit in m/s, then km/h.
- The DOPPLER frequency measurements and the calculation of the speed of displacement can be made for different speeds of the mobile unit.

# WAVEGUIDE PROPAGATION

## 4.10. DIRECTIVE COUPLER STUDY

### 4.10.1 Introduction

The manipulations that follow are designed to highlight the basic parameters of the couplers: coupling and directivity then the use of directive couplers to simply measure the standing wave ratio. These manipulations are valid both for the GCX 100 cross coupler and for the CDT 100 the directive hole coupler.

### 4.10.2 Theoretical reminders

#### ■ Directive coupler

To evaluate the microwave power conveyed by a transmission line connecting a transmitter to an antenna, we use a directive coupler combined with a detector and indicator assembly. This measurement consists in sampling part of the incident wave independently of the energy reflected by the transmission line load impedance.

The directive coupler consists of:

- a waveguide with the main arm in series in the transmission line.
- a waveguide section, the auxiliary arm adjacent to the main arm, used for sampling the coupled energy.

One end of the auxiliary arm is used for connecting a detector system and the other is closed on a matched load absorbing the residual energy transmitted, due to coupler non-infinite directivity.

These two waveguides are coupled by holes, slots, crosses, etc, depending on the type of coupler. The two essential characteristics of the coupler are coupling and directivity.

#### ■ Coupling

The main parameter of a directive coupler is coupling. Expressed in dB, it corresponds to the power ratio P1 picked up by the coupled output from the incident power P.

Coupling example :

Cross-type coupler coupling	> 20 dB to 60 dB
Hole coupler coupling	3 to 30 dB.

$$C(\text{dB}) = 10 \log \frac{P_1}{P}$$


fig : 1

Coupling can be 3, 10, 20, 30, 40, 60 dB, depending on the use.

There are low coupling value (in dB) couplers used for transmitting energy to RF mixer heads used in radar. These hybrid couplers halve the received energy and are called "dividers" or "splitters".

■ **Electromagnetic coupling**

In the case of coupling by holes in the axis of the long side of a waveguide, the distribution of the electric field in the auxiliary waveguide near the hole is similar to that which would be produced by an electric dipole whose moment is parallel to the electric field of the main waveguide and in the same direction.



fig : 2

equivalent dipole

Similarly, the magnetic field in an auxiliary waveguide behaves as if the hole were radiating like a magnetic dipole whose moment is parallel to the transverse magnetic field of the main waveguide, but this time, in the opposite direction.



fig : 3

equivalent dipole

It is the change of direction that cancels the wave in one of the directions of the coupler auxiliary arm, creating the coupler directivity.

■ **Directivity**

Directivity is also evaluated in dB. It expresses the relative quantity of coupled power P1 with respect to the return power P2 transmitted to the auxiliary arm. Therefore, it can be considered as a measurement of coupling quality. The higher the directivity, the better the quality of the coupler.

Directivity example:

- Cross-type coupler directivity ≈ 15 to 20 dB.
- Hole coupler directivity ≈ 20 to 40 dB.

$$D(dB) = 10 \log \frac{P1}{P2}$$


fig : 4

A coupler having high directivity is barely sensitive to the energy received by reflection on the load impedance of the transmission line.

### 4.10.3 Coupler experimentation

#### ■ Composition :

- 1 Cross coupler:	GCX 100	or 1 directive hole coupler:	CDT 100/20
- 1 Iris 20 and 30 :	IRIS 100	or 1 Iris 30 dB :	ICDT 100/30
- 1 oscillator :	OSG 100		
- 1 ferrite isolator :	ISO 100		
- 1 PIN modulator :	MOD 100		
- 1 micrometer attenuator:	ATM 100		
- 1 calibrated attenuator :	ATC 100		
- 1 waveguide/coaxial transition:	TGN 100		
- 1 measurement sensor:	DEN 100	or 1 thermocouple probe:	ST 613
- 2 matched loads:	CHG 100		
- 3 waveguide supports:	SUP 100		
- 1 measurement indicator. :	IR 205	or 1 powermeter:	MH 600
- 1 Gunn power supply:	CF 204		

#### ■ Starting up the bench

The bench consists of an OSG 100 Gunn diode oscillator. The microwave unit is modulated by 1 kHz square signals using a MOD 100 - PIN diode modulator. An ISO 100 ferrite isolator placed between the oscillator and the pin diode modulator prevents load variations caused by the latter from reacting on the oscillator. A set of 2 ATM 100 and ATC 100 attenuators controls the level and allows attenuation measurements to be made

- In succession, we will be studying the GCX 100 cross coupler and the CDT 100 hole coupler. These couplers are removable and an effort will be made to investigate the various forms of coupling openings by exchanging the coupling irises (IRIS 100 and ICDT 100/20 and /30).
- Assemble the microwave components forming the bench, connecting the IR 205 SWR indicator as indicated below.
- Power the CF 204 Gunn power supply and the IR 205 SWR indicator.
- Proceed with the measurements.

#### ■ Coupling measurement

Set up as indicated below (figure 5). (connect the coupler in compliance with the direction shown by the arrow).

The following components must be mounted on the coupler arm as mentioned below:

- 1) The TGN 100 waveguide-coaxial transition and the DEN 100 detector on the coupled output (2)<sup>1</sup>, [4]<sup>2</sup>.
- 2) The loads (3) and (4), [3].


fig : 5

<sup>1</sup> (x) : CGX100. Output

<sup>2</sup> [x] : CDT100. Output

x : see on the figure, the corresponding position.

- Set the ATC 100 attenuator to 0 dB and adjust the ATM 100 attenuators to obtain correct deflection on the galvanometer of IR 205.
- Read the indication given on IR 205.
- Then set up as indicated in the following figure (fig 6).

The following components must be mounted on the coupler arm as mentioned below:

- 1) The TGN 100 waveguide-coaxial transition and the DEN 100 detector on the coupled output (3)<sup>1</sup>, [3]<sup>2</sup>.
- 2) The loads (2) and (4), [4].


fig : 6

- Adjust the ATC 100 attenuator to bring the indication on the galvanometer of IR 205 to the mark obtained with the previous set up.
- On the dial of the ATC 100 attenuator, read the attenuation corresponding to the measured coupling value.

These manipulations can be carried out at several frequency points to check whether the coupling is constant as a function of the frequency.

<sup>1</sup> See notes, § 4.10.3 coupler expérimentation.

<sup>2</sup> See notes, § 4.10.3 coupler expérimentation.

■ **Directivity measurement**

Set up as shown below (fig 7):

The following components must be mounted on the coupler arm as mentioned below:

- 1) The TGN 100 waveguide-coaxial transition and the DEN 100 detector on the coupled output (4)<sup>1</sup>, [4]<sup>2</sup>.
- 2) The loads at (2) and (3) [1].


fig : 7

- Set the ATC 100 attenuator to 0 dB and read the indication on the galvanometer of IR 205.
- Locate the deflection obtained on the galvanometer.

Set up as indicated below (fig 8) :

The following components must be mounted on the coupler arm as mentioned below:

- 1) The TGN 100 waveguide-coaxial transition and the DEN 100 detector on the coupled output (3), [3].
- 2) The loads at (2) and (4), [4].


fig : 8

<sup>1</sup> See notes, § 4.10.3 coupler experimentation.

<sup>2</sup> See notes, § 4.10.3 coupler experimentation.

- Adjust the ATC 100 attenuator to bring the indication of the galvanometer to the mark obtained with the previous assembly and read directly the attenuation corresponding to the desired directivity.

These manipulations can be carried out at several frequency points to check whether the coupling is constant as a function of the frequency.

NOTE: It is necessary to use a sensitive measurement range on SWR indicator to obtain a legible deviation on the galvanometer .

## 4.11. STUDY OF FERRITE ISOLATOR AND CIRCULATOR

### 4.11.1 Introduction

The purpose of this series of manipulations is to highlight the installation and switching properties of the ferrite devices. They are described for the ferrite circulator.

Experimentation on a ferrite isolator is similar to that of a circulator but replacing in the setup the circulator by the isolator to be investigated.

### 4.11.2 Theoretical reminders

#### ■ General :

The ferrite consists of magnetic iron oxides to which are added other components such as Yttrium, Nickel or Manganese.

High resistivity and permeability allow microwave use.

The magnetic properties of the ferrites are linked directly to the magnetic moments of the spins themselves due to the fact that each atom includes unmatched electrons.

If the ferrite is placed in a magnetic field, the orientation of the magnetic moments of the atoms is modified.

The magnetization created in this way is proportional to the outside magnetic field until saturation is achieved.

The permeability of a saturated ferrite changes according to the magnetization direction: this property is used in microwave applications to obtain components having directional characteristics, like circulators or isolators.

#### ■ The circulator :

The circulator consists of three waveguide sections connected to form a Y and terminating in connecting flanges at their free ends.

The bar of ferrite and the middle of the junction is used for directing the microwave energy in a privileged direction, when magnetized by an outside magnetic field "H" : see the figure below


fig : 1

Accordingly, microwave energy injected at the input (1) is encountered at the output (2) only.

Similarly, by circular changeover, the energy injected at (3) or (2) is found at (1) or (3).

For instance, the circulator isolates the receiver from the transmitter when it is coupled with the same antenna.

#### ■ The isolator :

The ferrite isolator consists of a straight waveguide section both ends of which are provided with a connecting flange.

A (cylindrical) piece of ferrite, placed on the inside, is magnetized by an outside magnetic field "H"..

The ferrite isolator allows the microwave energy to pass in one direction but opposes it, by absorbing it, in the other direction.

A circulator can be used as an isolator when it is necessary to transmit microwave energy from input (1) to output (2), as long as the output (3) is closed on a matched load.

The ISO 100 isolator supplied with the BDH R 100 bench is of this type.

### 4.11.3 Experimentation

#### ■ Composition :

- |                                    |         |        |                              |
|------------------------------------|---------|--------|------------------------------|
| - 1 Gunn oscillator :              | OSG 100 |        |                              |
| - 1 ferrite isolator :             | ISO 100 |        |                              |
| - 1 PIN modulator :                | MOD 100 |        |                              |
| - 1 micrometer attenuator :        | ATM 100 |        |                              |
| - 1 calibrated attenuator:         | ATC 100 |        |                              |
| - 1 ferrite circulator :           | CIR 100 | and/or | 1 isolator: ISO 100          |
| - 3 fastening adapters :           | AFR 100 |        |                              |
| - 1 matched load :                 | CHG 100 |        |                              |
| - 1 waveguide/coaxial transition : | TGN 100 |        |                              |
| - 1 measurement sensor :           | DEN 100 | or     | 1 thermocouple probe: ST 613 |
| - 3 waveguide supports :           | SUP 100 |        |                              |
| - 1 measurement indicator :        | IR 205  | or     | 1 powermeter: MH 600         |
| - 1 Gunn power supply:             | CF 204  |        |                              |

#### ■ Starting up the bench

The components used for performing the experiments that follow are, in their great majority, BDH R100 basic bench components.

- Set up the microwave components and connect the CF 204 power supply, the IR 205 SWR indicator, or the MH 600 + ST 613 powermeter and ST 613 thermocouple probe as shown below (fig. 2).
- Power up the CF 204 Gunn power supply, IR 205 SWR indicator or the powermeter.
- Proceed with the measurements.

#### ■ Insertion loss measurement

In the following case, use either a CIR 100 circulator (or an ISO 100 isolator).

Set up as shown in figure 2.

(the unused output (3) of the circulator is closed on a matched load CHG 100).


fig : 2

- Power up the OSG 100 Gunn oscillator and adjust the operating frequency to approximately 9 GHz.
- Using the variable attenuator ATM 100, adjust the power measured by the MH 600 powermeter to obtain 0 dB (relative mode)
- Remove the CIR 100 circulator and measure the power available at the output of the variable attenuator ATM 100.
- Read directly in dB on MH 600 the circulator insertion loss.
- Transfer this value to a graph representing the insertion loss of the circulator according to the frequency and make the measurement again for other frequencies included between 8.5 GHz and 9.6 GHz.
- Begin the experiment in the direction (3) --> (1), then (2) --> (3).

■ **Insulation measurement**  
Set up as shown in figure 3.


fig : 3

- Set the precision attenuator ATC 100 to 0dB ;
- Power up the OSG 100 Gunn oscillator and adjust the operating frequency to approximately 9 GHz.
- Tune the 1 kHz modulation frequency to obtain the best possible measurement sensitivity on the IR 205 SWR indicator or MH 600 powermeter.
- Adjust the attenuation generated by the variable attenuator ATM 100 to obtain a reference indication  $I_0$  on IR 205 (or MH 600).
- Remove the circulator and adjust the attenuator ATC 100 to obtain the deflection  $I_0$  again.
- The attenuation value given by ATC 100 corresponds to the isolation of the circulator CIR 100.
- Transfer this value to a graph representing the insertion loss of the circulator according to the frequency and make the measurement again for other frequencies included between 8.5 GHz and 9,6GHz.
- Begin the experiment in the direction (3) --> (2), then (1) --> (3).

## 4.12. DIELECTRIC STUDY

### 4.12.1 Introduction

The order of magnitude of the dielectric constant of most materials habitually used for microwave applications is known but it is sometimes necessary to measure this constant accurately. In addition, the dielectric constant of some materials, like composites or ceramics, is not necessarily known and has to be estimated by experiment.

The same applies to the loss tangent.

### 4.12.2 Theoretical reminders

#### ■ Dielectric constant of the dielectric medium

The electric fields  $\vec{E}$  and magnetic fields  $\vec{H}$  propagate in a perfect dielectric (conductivity  $\sigma = 0$ ) that is not limited (free space), in conformity with the following equations, deduced from those of MAXWELL.

$$\Delta \vec{E} - \mu\epsilon \frac{\delta^2 \vec{E}}{\delta t^2} = 0$$

$$\Delta \vec{H} - \mu\epsilon \frac{\delta^2 \vec{H}}{\delta t^2} = 0$$

These equations, referred to as wave equations, are obtained by applying conventional mathematical identity:

$$\text{rot} \cdot \text{rot} \cdot \vec{A} = \text{grad} \cdot \text{div} \cdot \vec{A} - \Delta \vec{A}$$

with the following system of equations:

$$\text{rot} \cdot \vec{E} = -\mu \frac{\delta \vec{H}}{\delta t}$$

$$\text{rot} \cdot \vec{H} = -\epsilon \frac{\delta \vec{E}}{\delta t}$$

where  $\text{div} \cdot \vec{E} = 0$  and  $\text{div} \cdot \vec{H} = 0$

MAXWELL's equation in the absence of free loads and conduction current, in a non conductive medium.

In the harmonic configuration the time dependence of the fields results in a factor  $e^{j\omega t}$ ; these equations are written:

$$\begin{aligned} \Delta \vec{E} + k^2 \vec{E} &= 0 \\ \Delta \vec{H} + k^2 \vec{H} &= 0 \end{aligned}$$

with  $k^2 = \omega^2 \epsilon\mu$

#### ■ Loss angle

For a conductivity dielectric that is not zero,  $\sigma$ , the MAXWELL-AMPERE equation in the harmonic configuration is written:

$$\text{rot} \cdot \vec{H} = (\sigma + j\omega\epsilon) \vec{E}$$

hence :

$$\sigma \ll \omega\epsilon$$

The second member is written :

$$j\omega \left( \epsilon - j \frac{\sigma}{\omega} \right) \vec{E}$$

If we set:

$$\tilde{\varepsilon} = \varepsilon' - j\varepsilon'' = \varepsilon - j \frac{\sigma}{\omega}$$

$\tilde{\varepsilon}$  is a complex magnitude giving the form :

$$\text{rot. } \vec{H} = j \omega \tilde{\varepsilon} \vec{E}$$

The  $\text{tg } \delta = \varepsilon'' / \varepsilon'$  ratio is the loss factor or tangent of loss angle  $\delta$ .

In current dielectrics, this loss factor is on the order of:  $10^{-3}$  to  $10^{-4}$

### 4.12.3 Experimenting with dielectrics

#### ■ Composition :

- 1 Gunn oscillator :	OSG 100		
- 1 ferrite isolator :	ISO 100		
- 1 PIN modulator :	MOD 100		
- 1 straight waveguide :	GD 100/180		
- 1 slotted line :	LAF 100		
- 1 matched load :	CHG 100		
- 1 Wavemeter :	OND 100		
- 1 dielectric kit :	KED 100		
- 1 micrometer attenuator :	ATM 100		
- 1 calibrated attenuator :	ATC 100		
- 3 waveguide supports :	SUP 100		
- 1 measurement sensor :	DEN 100	or 1 thermocouple probe:	ST 613
- 1 Gunn power supply :	CF 204		
- 1 measurement indicator. :	IR 205	or 1 powermeter:	MH 600

#### ■ Dielectric study bench activation

An OSC 100 oscillator generates the microwave. It is modulated (1 kHz) by a modulator MOD 100. An ISO 100 ferrite isolator placed between the oscillator and the modulator prevents the load variations caused by the latter from reacting on the oscillator. An ATM 100 attenuator checks the level and an OND 100 wavemeter allows the frequency to be measured.

LAF 100 slotted line is used for determining the SWR and the electric field maximum and minimum positions in the waveguide.

A straight section of waveguide, GD 100/180, terminating by a CC 100 short-circuit is used for inserting the samples to be tested.

#### 4.12.3.1. Measurement of dielectric constant (small parallelepipedic samples)

Method of short-circuit line (short parallelepipedic samples)

Set up the microwave components and connect the slotted line detector LAF 100 to the IR 205 indicator.(fig. 1)

Power the CF 204 Gunn power supply and the IR 205 SWR indicator.

Proceed to measurement.


fig : 1

The sample is a parallelepipedic sized to match the straight section of the waveguide. Adjustment is made particularly carefully along the longest side of the waveguide. Finally, the sample must be brought into perfect contact with the short-circuit terminating the waveguide

#### Measurement method

##### CAUTION !

Prevent the complete depression of the slotted line plunger (leave a 4 mm clearance) so as to prevent a high SWR from being combined with the phenomenon to be studied.

With the sample absent, locate accurately the 0 position of the standing wave. The introduction of a sample will cause this point to move in the direction of the short-circuit. Let  $\Delta$  be the measurement movement value.

The parameter U is obtained by means of the following equation:

$$U = \frac{\Delta + d}{\lambda_g} \quad (1)$$

where d is the thickness of the sample and  $\lambda_g$  the wavelength in the waveguide in the absence of a dielectric.

Then V is given by the equation:

$$\frac{\text{tg } 2\pi V}{V} = \frac{\lambda_g}{d} \text{tg } 2\pi U \quad (2)$$

and K, the dielectric constant is given by :

$$V = \frac{d\sqrt{K - p}}{\lambda} \quad (3)$$

where

$$p = \left[ \frac{\lambda_0}{\lambda_c} \right]^2$$

with:  $\lambda_0$  = the wavelength in the air

$\lambda_c$  = the waveguide cut-off wave length without a dielectric

The equation (2) is expressed graphically by a curve:  $\frac{\text{tg } x}{x}$

##### Measurement accuracy

The quantities d and  $\lambda$  can be obtained accurately without excessive difficulties using a calibre gauge (for d) and a wavemeter (for  $\lambda$ ).

The d/ $\lambda$  ratio is therefore known. The error U depends on the precision of the measurement of  $\Delta$  and  $\lambda_g$

#### 4.12.3.2. Slotted live method (slotted samples) parallelepipedic

The composition of the bench is the same as previously (fig. : 1)

#### Measurement method

##### CAUTION !

Prevent the complete depression of the slotted line plunger (leave a 4 mm clearance) so as to prevent a high SWR from being combined with the phenomenon to be studied.

Using the LAF 100 slotted line, locate in the position of the maxima so as to deduce  $\lambda$  in the air then insert the sample into the slotted line and make new measurements so as to obtain  $\lambda_{gd}$  in the dielectric.

It is impossible to calculate :  $\varepsilon = \left[ \frac{\lambda_g}{\lambda_{gd}} \right]^2$

With:  $\lambda_{gd}$  = wavelength with dielectric

#### 4.12.3.3. Measurement of loss angle (long beveled sample)

Modify the previous setup (fig: 1) to make it compliance with the setup shown below (fig. :2)


fig :2

#### Measurement method

The sample is inserted into the waveguide whose perception is a parallelogram enabling the wave to enter or exit the dielectric by gradual impedance variation and generating only very slight reflection. It is essential that the sample bevels on between the two larger surfaces. Otherwise, propagation in mode  $TE_{02}$  would exist in the waveguide filled with the dielectric.

After inserting the sample, the wave is attenuated in the ratio A and B is measured directly in dB.

The value is measured by means of the calibrated attenuator ATC 100, making it possible to return to the level read on the powermeter before the insertion of the sample.

We obtain  $tg\delta$  through the equation :

$$tg \delta = \frac{\lambda \sqrt{K - p \log_e A}}{ndk}$$

in which d is the length of the sample (equivalent to that of a parallelepipedic section having the same volume).

This equation only applies to low values of  $tg\delta$ .

The value A is the ratio between the electric field before and after passing through the sample measured in dB.

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